



I am your outdoor posterior walking aid.

Just like my brother, I am always with you when you really want to get moving, and will offer you the best-possible support to do so. Thanks to my outdoor outfit, I am happiest outdoors, and am the perfect companion for a walk: my large wheels will roll along any surface and my drum and parking brakes provide additional security.



Optional anti-tip mechanism
protects you from tipping back.

Frame colours size 1



Frame colours size 2



Frame colours size 3





	Size 1	Size 2	Size 3
Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm
Height of forearm rests	46 - 64 cm	63 - 87 cm	79 - 103 cm
Handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm
Seat height	33 cm	45 cm	60 cm
Total width	56 cm	66 cm	73 cm
Total length	66 cm	75 cm	82 cm
Frame height (folded)	25 cm	32 cm	34 cm
Transport dimensions (L x W x H)	66 x 56 x 25 cm	75 x 66 x 32 cm	82 x 73 x 34 cm
Max. load	30 kg	50 kg	70 kg
Weight	7.1 kg	8.3 kg	9.8 kg

Did you know that I am one of the few walking trainers which has an outdoor approval? With larger and wider wheels and a drum and parking brake, I am well-equipped for big adventures outdoors. Brilliant!



You can find more information on accessories in the price list.

