

GMFCS
Level III-IV

kg
max. 100



46–124 cm

29–113 cm



I am your versatile posterior walking aid.

I am always there for you – even from your first year of life. Walking is my biggest hobby. It does not matter where you want to go; I will follow you. If you support yourself on my handles, you will quickly notice how I take on your body weight. In addition, my special abducted frame shape supports you so that you do not swing outwards when walking. My pelvic support pad ensures that you walk upright and straight. Cool, aren't I?!





Combined pelvic and buttock pad
can become a backrest in one single easy step.

Forearm rests
with handles.

Padded seat
can easily be folded back.

Aluminium frame
can easily be folded up.

Direction locks
facilitate
walking straight.

Backwards stop
prevents inadvertent
rolling back.

EVA guide wheels
make the walking aid manoeuvrable.



Abducted frame shape
supports you when walking
and prevents you from
swinging outwards.



Aluminium frame
can easily be folded up via
folding joints.

Frame colour
size 0 mini

● Lazurite blue
metallic

Frame colours
size 1

● Telemagenta

● Sky blue

Frame colours
size 2

● Lazurite blue
metallic

● Traffic red

Frame colours
size 3 and size 4

● Anthracite
metallic

○ Pure white



	Size 0 mini*	Size 1	Size 2	Size 3	Size 4
Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm
Height of forearm rests	-	46 - 70 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm
Handle spacing	22 - 36 cm		30 - 44 cm	38 - 52 cm	43 - 57 cm
Seat height	-	33 cm	45 cm	60 cm	70 cm
Total width	52 cm	51 cm	61 cm	68 cm	75 cm
Total length	58 cm		70 cm	77 cm	87 cm
Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm
Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm
Max. load	20 kg	30 kg	50 kg	75 kg	100 kg
Weight	3.9 kg	7.1 kg	8.3 kg	8.5 kg	10.2 kg

* Frame cannot be folded up



You can find more information on accessories
in the price list.

