



Pricelist 2021

**Hello!**

we are Miriam and Torsten Schuchmann. Our family-run company is now in its second generation, our company has over 30 years of experience in producing and selling childrens' rehabilitation equipment.

It is a matter close to our hearts to support children with individual requirements as optimally as possible – using products which are just as individual as you are.

What inspires us.

The reason we have focused our efforts on childrens' rehabilitation equipment has its origins in our own family. We have first-hand experience and completely understand the challenges associated with a severely disabled person with multiple disabilities.

When Miriam's sister was young the development of rehabilitation equipment was at a very early stage. Her parents knew what they could improve to facilitate her participation in everyday life, and a business concept was born.

What we inspire.

Our aim is to support children with disabilities so that they can become more mobile and independent. At the same time, our task has more far-reaching implications. Meanwhile our 150 man team works closely with specialist dealers, doctors and therapists so that we can truly understand the requirements of the children and their families. An extremely moving process for all those involved.

Our products represent more than just products. They become companions for life. And last of all, we want to get the market moving – with real innovations in rehabilitation equipment which fulfil the maximum demands on functionality, design and safety.

If you have any questions or suggestions on general topics or on our products, we – and of course the entire Schuchmann Team – are happy to help at any time.

Seating.



smilla.
04–07



madita-fun.
08–13



madita-fun mini.
14–16



madita hd.
18–21



tina.
22–24

Transport.



new

carrot 3.
26–29



malte.
30–32



malte outdoor.
34–35



marcy.
36–37

Walking.

Standing.



till.
40–45



todd.
46–51



tim.
52–55

Hygiene.



boris.
58–59



basti.
60–61



bibi.
62

General guidelines.



**Cleaning +
disinfection.**
64–65

Contents.



smilla.
Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____


Reference:
The seating aid for everyday use.

- Sturdy wooden frame with integrated aluminium profiles
- Incl. height adjustable backrest and adjustable back angle
- Incl. seat and back pad
- Incl. height adjustable armrests
- Incl. height, angle and depth adjustable seat and footplate
- Incl. positioning castors and pelvic harness
- Material: varnished beech


Dimensions - basic model

		Size 0 Nursery	Size 0	Size 1 Nursery	Size 1	Size 2
A	Seat width	17 - 28 cm	17 - 28 cm	20 - 35 cm	20 - 35 cm	30 - 40 cm
B	Seat depth	16 - 26 cm	16 - 26 cm	25 - 35 cm	25 - 35 cm	30 - 40 cm
C	Max. seat height	45 cm	65 cm	45 cm	65 cm	65 cm
D	Back height	25 cm	25 cm	29 cm	29 cm	34 cm
E	Back angle	(-) 3° - 17,5°				
F	Tilt-in-space	(-) 5°, 0°, 7,5°, 15°				
G	Foot angle	(-) 5°, 0°, 7,5°, 15°				
H*	Min. footplate height	Standard	13 cm			
		with 50 mm dual castors	18 cm			
		with 75 mm dual castors	21 cm			
H*	Min. footplate height (with quick angle adjustment)	Standard	8 cm			
		with 50 mm dual castors	13 cm			
		with 75 mm dual castors	16 cm			
	Max. load	30 kg	30 kg	45 kg	45 kg	60 kg
	Weight	9.2 kg	9.7 kg	9.6 kg	10.2 kg	12 kg

* Max. adjustable lower leg length = selected / required seat height - Min. footplate height (H)





Basic model		
smilla. size 0	smilla. size 1	smilla. size 2
Price	Price	Price
☐ Size 0 Nursery (13 00 100/UK) £ 620	☐ Size 1 Nursery (13 01 100/UK) £ 620	☐ Size 2 (13 02 000/UK) £ 720
☐ Size 0 (13 00 000/UK) £ 620	☐ Size 1 (13 01 000/UK) £ 620	

Frame colours	
Beech natural	Beech white
	
Price	Price
☐ 92 00 300 incl.	☐ 92 00 301 £ 55

Covers					
Pink (No. 115)	Azur (No. 110)	Green (No. 89)	Dark grey (No. 92)	Flexi Blue (No. 44)	Name embroidery on the back pad please fill out the form on page 25! Embroidery of flexi fabric is not available.
					
Price	Price	Price	Price	Upcharge	
☐ 91 00 115 incl.	☐ 91 00 110 incl.	☐ 91 00 089 incl.	☐ 91 00 092 incl.	☐ 91 00 044 £ 145	

Accessories						
Headrest – padded and height adjustable		Lateral pads – padded and adjustable		Thigh guides – padded and adjustable		
						
Adjustable height (in 3 cm steps)					Length	Height
for size 0+1	12 – 21 cm	for size 0	13 – 19 cm	Width 1	17.5 cm	7 cm
for size 2	20.5 – 32.5 cm	for size 1	20 – 26 cm	Width 2	20 cm	8.5 cm
		for size 2	24 – 30 cm		22.5 cm	10 cm
Price		Price / pair		Price / pair		
☐ for size 0 (13 00 001)	£ 160	☐ Width 1 (13 01 002)	£ 100	☐ Size 0 (13 00 004)	£ 135	
☐ for size 1 (13 01 001)	£ 160	☐ Width 2 (13 01 003)	£ 100	☐ Size 1 (13 01 004)	£ 135	
☐ for size 2 (13 02 001)	£ 160			☐ Size 2 (13 02 004)	£ 135	

Configuration possibilities - seat	
Standard seat pad – incl. cover (standard configuration, if nothing else stated)	Seat pad – incl. cover with anatomic leg guidance
	
Price	Upcharge
☐ for size 0 (13 00 021)	☐ for size 0 (13 00 022) £ 75
☐ for size 1 (13 01 021)	☐ for size 1 (13 01 022) £ 75
☐ for size 2 (13 02 021)	☐ for size 2 (13 02 022) £ 75

Accessories

Quick angle adjustment of the footplate and seat – allows quick and easy angle adjustment of the footplate and seat to (-) 5°, 0°, 7.5° and 15° by straps (not possible in combination with quick height adjustment)



Upcharge / set

- ☐ for size 0 **(13 00 090)**
- ☐ for size 1 **(13 01 090)**
- ☐ for size 2 **(13 02 090)**

£ 165
£ 165
£ 165

Quick height adjustment for footplate, seat and back – allows quick and easy height adjustment of the footplate, seat and back (not possible in combination with quick angle adjustment)



Upcharge / set

- ☐ **13 00 020**

£ 42

Accessories

Wooden therapy tray – plug-on



	Width	Depth
for size 0	48 cm	37 cm
for size 1	52 cm	39 cm
for size 2	57 cm	39 cm

Price

- ☐ for size 0 **(19 00 005)**
- ☐ for size 1 **(19 01 005)**
- ☐ for size 2 **(19 02 005)**

£ 140
£ 140
£ 140

Large wooden therapy tray – slide on, with edge



	Width	Depth
for size 0	48 cm	45 cm
for size 1	52 cm	49 cm
for size 2	57 cm	49 cm

Price

- ☐ for size 0 **(19 00 009)**
- ☐ for size 1 **(19 01 009)**
- ☐ for size 2 **(19 02 009)**

£ 165
£ 165
£ 165

Abduction wedge – padded and quick release



Price

- ☐ for size 0 **(19 00 008)**
- ☐ for size 1+2 **(19 01 008)**

£ 98
£ 98

Accessories

Dual castors – 4 pieces



Price / set

- ☐ 50 mm – 2 with locking devices **(13 01 006)** £ 160
- ☐ 50 mm – all with locking devices **(13 01 016)** £ 185
- ☐ 75 mm – 2 with locking devices **(13 00 006)** £ 160
- ☐ 75 mm – all with locking devices **(13 00 016)** £ 185

Harness

The **smilla.** harness is assembled with the following harness assembly sets.
If a colour is not selected, we will use black as the standard colour for your order.

Harness assembly kit for the seat
consisting of 2 brackets with clip buckle and 2 clip buckles (drilled)



Price / set

- ☐ 20 mm **(19 00 010)**
- ☐ 25 mm **(19 01 010)**

incl.
incl.

Harness assembly kit for the back
consisting of belt holder with T-groove incl. 2 clip buckles and 2 clip buckles (drilled)



Price / set


- ☐ 20 mm **(19 00 011)**
- ☐ 25 mm **(19 01 011)**

£ 55
£ 55

Harness


Chest harness – made of elastic, breathable material, three straps with clips, sewn strap to the left shoulder

	Size 0	Size 1	Size 2	Size 3
A	13cm	13cm	18cm	18cm
B	17.5cm	19cm	20cm	21cm
C	15cm	15cm	19cm	28cm
D	3.5cm	5cm	5cm	5cm
E	25.5cm	28cm	34cm	44cm
Strap width	2cm	2cm	2.5cm	2.5cm
Price Colour				
☐ Size 0 (1300025) ☐ Size 1 (1301025) ☐ Size 2 (1302025) ☐ Size 3 (1303025) £ 99 £ 99 £ 99 £ 99 ☐ Black (9100040) ☐ Light grey (9100049)				



Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips


	Size 0	Size 1	Size 2
A	11cm	11cm	18cm
B	21cm	21cm	22cm
C	12cm	12cm	15cm
D	2.5cm	3.5cm	3.5cm
E	25.5cm	28cm	34cm
Strap width	2cm	2cm	2.5cm
Price Colour			
☐ Size 0 (1300026) ☐ Size 1 (1301026) ☐ Size 2 (1302026) £ 119 £ 119 £ 119 ☐ Black (9100040) ☐ Light grey (9100049)			



Harness


T-shape groin harness – made of elastic, breathable material, with two clips at front, straps sewn at back

	Size 0	Size 1	Size 2
A	23cm	25cm	28cm
B	15cm	17cm	20cm
C	33cm	35cm	38cm
D	24cm	26cm	28cm
Strap width	2cm	2cm	2.5cm
Price Colour			
☐ Size 0 (1300027) ☐ Size 1 (1301027) ☐ Size 2 (1302027) £ 88 £ 88 £ 88 ☐ Black (9100040) ☐ Light grey (9100049)			



Groin harness (standard) – made of elastic, breathable material, with two clips at front, straps sewn at back

	Size 1	Size 2	Size 3
A	24cm	29cm	33cm
B	14cm	17cm	22cm
C	35cm	43cm	49cm
Strap width	2cm	2.5cm	2.5cm
Price Colour			
☐ Size 1 (1301028) ☐ Size 2 (1302028) ☐ Size 3 (1303028) £ 88 £ 88 £ 88 ☐ Black (9100040) ☐ Light grey (9100049)			



Harness

Pelvic harness – with pad and adjustable clip on both sides (standard configuration = pelvic harness size 0 for **smilla**. size 0, pelvic harness size 1 for **smilla**. size 1 and pelvic harness size 2 for **smilla**. size 2, if nothing else stated)


	Size 0	Size 1	Size 2
Pad width	4cm	4cm	5cm
Pad length per side	16cm	18cm	28cm
Strap length per side	30cm	30cm	40cm
Strap width	2cm	2.5cm	2.5cm
Price Colour			
☐ Size 0 (1300029) ☐ Size 1 (1301029) ☐ Size 2 (1302029) incl. incl. incl. ☐ Black (9100040) ☐ Light grey (9100049)			



Harness

4-point pelvic harness – with pad and adjustable clip on both sides

	Size 0	Size 1	Size 2
Pad width	4cm	4cm	5cm
Pad length per side	16cm	18cm	28cm
Strap length per side	30cm	30cm	40cm
Strap width	2cm	2.5cm	2.5cm
Upcharge Colour			
☐ Size 0 (1300030) ☐ Size 1 (1301030) ☐ Size 2 (1302030) £ 32 £ 32 £ 32 ☐ Black (9100040) ☐ Light grey (9100049)			



madita-fun.

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____

Reference:
The modern and modular therapy chair.

- Height adjustment via gas pressure spring 1200 N
- Adjustable pelvic guidance
- Incl. push handle, height, angle and width adjustable armrests
- Incl. 4-point pelvic harness, assembly kits for the seat and back
- Back unit in short / long standard design
- Incl. folding footplate, knee angle, foot angle and height adjustable
- 4 castors (75 mm), 2 with locking devices
- Adjustment options:
 - tilt-in-space
 - back angle
 - seat depth
 - seat width



Basic model

madita-fun. size 0	madita-fun. size 1	madita-fun. size 1b	madita-fun. size 2
Price	Price	Price	Price
☐ 74 00 100/UK £ 1,595	☐ 74 01 000/UK £ 1,595	☐ 74 01 100/UK £ 1,660	☐ 74 02 000/UK £ 1,760

Frame colours

Anthracite metallic for size 1, 1b & 2		Telegrey 4 (white grey) for size 0, 1, 1b & 2	
Price	Price	Price	Price
☐ 92 00 014 incl.	☐ 92 00 011 incl.	☐ 92 00 011 incl.	☐ 92 00 011 incl.

Presettings



Presettings on the therapy chair at Schuchmann.
Please provide us with the following body measurements:

A Axilla height	_____	D Pelvic width	_____
B Chest width	_____	E Seat depth	_____
C Back height	_____	F Lower leg length	_____


☐ **74 00 080**

Price

incl.

Dimensions - basic model

		Size 0	Size 1	Size 1 b	Size 2
A	Seat depth*	15 - 25 cm	21 - 30 cm	26 - 36 cm	30 - 42 cm
B	Seat width	16 - 26 cm	20 - 30 cm	20 - 30 cm	28 - 38 cm
C	Back height with short back	32 - 37 cm	35 - 43 cm	42 - 52 cm	42 - 52 cm
C	Back height with long back	35 - 43 cm	42 - 52 cm	53 - 63 cm	53 - 63 cm
D	Seat height	25 - 62 cm	25 - 62 cm	25 - 62 cm	25 - 62 cm
	Seat height with electric height adjustment	32 - 62 cm	32 - 62 cm	32 - 62 cm	32 - 62 cm
	Seat height with hydraulic height adjustment	27 - 63 cm	27 - 63 cm	27 - 63 cm	27 - 63 cm
E	Lower leg length	17 - 25 cm	20 - 32 cm	28 - 40 cm	26 - 49 cm
F	Total width	54 cm	54 cm	64 cm	64 cm
G	Length of base frame	61 cm	61 cm	75 cm	75 cm
	Total length with combi subframe	76 cm	76 cm	86 cm	86 cm
	Tilt-in-space	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°	(-)8° - 37°
	Back angle	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°	(-)5° - 25°
	Max. load	30 kg	30 kg	45 kg	50 kg
	Weight	19 kg	20 kg	22 kg	24 kg



* reduced by 1 cm in combination with adjustable pelvic guide

Covers

Dark grey/ Pink-Telemagenta (No. 92 + 114) 	Dark grey/ Azur (No. 92 + 110) 	Dark grey/ Deep orange (No. 92 + 107) 	Dark grey/ Light blue (No. 92 + 106) 
Price	Price	Price	Price
☐ 91 00 222	☐ 91 00 225	☐ 91 00 213	☐ 91 00 214
incl.	incl.	incl.	incl.

Covers

Black/ Dark grey (No. 105 + 92) 	Jeans (No. 116) 
Price	Price
☐ 91 00 205	☐ 91 00 116
incl.	incl.

Name embroidery

Name embroidery on the back pad
please fill out the form on page 25!




Embroidery of flexi fabric is not available.

Coated covers

Flexi black/ Light grey (No. 43 + 50) 	Flexi black/ Red (No. 43 + 45) 	Flexi blue (No. 44) 
Upcharge	Upcharge	Upcharge
☐ 91 00 226	☐ 91 00 227	☐ 91 00 044
£ 152	£ 152	£ 152


Configuration options - back (Standard design or with integrated upper arm guide)

Back – incl. pad and cover	for size 0	for size 1	for size 1b	for size 2
Short Back: Long Back:	35 – 43 cm 42 – 52 cm	35 – 43 cm 42 – 52 cm	42 – 52 cm 53 – 63 cm	42 – 52 cm 53 – 63 cm
	Upcharge / piece	Upcharge / piece	Upcharge / piece	Upcharge / piece
Short standard design: Long standard design:	☐ 7300 102 incl. ☐ 7300 103 incl.	☐ 7301 202 incl. ☐ 7301 003 incl.	☐ 7301 102 incl. ☐ 7301 103 incl.	☐ 7302 002 incl. ☐ 7302 003 incl.
Short with upper arm guide: Long with upper arm guide:	☐ 7300 129 £ 47 ☐ 7300 030 £ 47	☐ 7301 229 £ 47 ☐ 7301 030 £ 47	☐ 7301 129 £ 47 ☐ 7301 130 £ 47	☐ 7302 029 £ 53 ☐ 7302 030 £ 53



Configuration options - back

Back unit – fixed (standard configuration = with push handle, if nothing else stated)	Back unit – removable	Pre-adjustable limitation of the max. backrest angle to determine the individual maximum back angle
		
Price	Price	Price / Set
☐ for size 1 with push handle (7301004) incl. ☐ for size 2 with push handle (7302004) incl. ☐ for size 0+1 with loop (7301053) incl.	☐ for size 1 with push handle (7301006) £ 242 ☐ for size 2 with push handle (7302006) £ 242	☐ 7301 111 £ 39

Headrests

Headrest (standard) padded, incl. cover, height and angle adjustable	Headrest (shell-shaped) padded, incl. cover, height and angle adjustable	Headrest (lateral guide) padded, incl. cover, height and angle adjustable	Headrest (lateral guide, angle adjustable) padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable																
																			
Price	Price	Price																	
☐ for size 0 (7300 108) £ 137 ☐ for size 1 (7301 008) £ 137 ☐ for size 2 (7302 008) £ 137	☐ for size 0+1 (7301 110) £ 195 ☐ for size 2 (7302 110) £ 195	☐ for size 0+1 (7301 009) £ 185 ☐ for size 2 (7302 009) £ 185	<table border="1"> <thead> <tr> <th></th><th>A</th><th>B</th><th>C</th></tr> </thead> <tbody> <tr> <td>Size 1</td><td>13.5 cm</td><td>10.5 cm</td><td>10 cm</td></tr> <tr> <td>Size 2</td><td>15 cm</td><td>10.5 cm</td><td>13.5 cm</td></tr> <tr> <td>Size 3</td><td>17 cm</td><td>12.5 cm</td><td>16.5 cm</td></tr> </tbody> </table> ☐ Size 1 (7301063) £ 348 ☐ Size 2 (7302063) £ 348 ☐ Size 3 (7303063) £ 348		A	B	C	Size 1	13.5 cm	10.5 cm	10 cm	Size 2	15 cm	10.5 cm	13.5 cm	Size 3	17 cm	12.5 cm	16.5 cm
	A	B	C																
Size 1	13.5 cm	10.5 cm	10 cm																
Size 2	15 cm	10.5 cm	13.5 cm																
Size 3	17 cm	12.5 cm	16.5 cm																

new



A = Width contact surface
B = Height headrest
C = Width lateral guide

Accessories

Universal adapter for tailored headrests – 15x15mm	Mounting point – similar Whitmyer bracket M2100-WB	Trigger lever – mounted directly on the gas spring	Trigger lever – mounted on the push bar/handle (standard configuration, if nothing else stated – not possible with removable back unit)
7301031 for sizes 0, 1, 1b and 2 (with short back) 7302031 for size 2 (with long back)			
Price		Price	Price
☐ 7301031 £ 40 ☐ 7302031 £ 40	☐ 7303031 £ 49	☐ 7400085 incl.	☐ 7401085 incl.

Tilt-in-space adjustment

Configuration options - lateral pads (Step 1: Select the bracket)				
Chest width – mounted on madita-fun. size 0	Chest width – mounted on madita-fun. size 1 + 1b	Chest width – mounted on madita-fun. size 2	Bracket	Upcharge / pair
13 – 20 cm	16 – 23 cm	26 – 32 cm	☐ Fixed (7300012) ☐ Flip away (7300013)	incl. £ 137
20 – 27 cm	23 – 30 cm	32 – 38 cm	☐ Fixed (7301012) ☐ Flip away (7301013)	incl. £ 137
___ – ___ cm	___ – ___ cm	___ – ___ cm	☐ Fixed and tailor made (7309012) ☐ Flip away and tailor made (7309013)	£ 47 £ 188

(Step 2: Select the lateral pads)

	Lateral pads size 0 padded, incl. cover (standard configuration for madita-fun. size 0, if nothing else stated) Height: 7 cm Depth: 9 cm	Lateral pads size 1 padded, incl. cover (standard configuration for madita-fun. size 1, if nothing else stated) Height: 8 cm Depth: 13 cm	Lateral pads size 2 padded, incl. cover (standard configuration for madita-fun. size 2, if nothing else stated) Height: 11 cm Depth: 16 cm	Tailored lateral pads padded, incl. cover Height: ____ cm Depth: ____ cm
	Price / pair	Price / pair	Price / pair	Price / pair
	☐ 73 00 011 £ 58	☐ 73 01 011 £ 58	☐ 73 02 011 £ 58	☐ 73 09 011 £ 115

Thigh guides

	Thigh guides for size 0 padded, incl. cover, adjustable angle, depth and width 10 cm = Standard configuration for size 0	Thigh guides for size 1 + 1 b padded, incl. cover, adjustable angle, depth and width 12 cm = Standard configuration for size 1 + 1 b	Thigh guides for size 2 padded, incl. cover, adjustable angle, depth and width 17 cm = Standard configuration for size 2
	Price / pair	Price / pair	Price / pair
	☐ 10 cm (73 00 114) £ 100 ☐ 12 cm (73 00 115) £ 100	☐ 12 cm (73 01 114) £ 100 ☐ 17 cm (73 01 115) £ 100	☐ 17 cm (73 02 114) £ 100 ☐ 22 cm (73 02 115) £ 100

Thigh guides with knee support

Thigh guides with knee-support for size 0 10 cm long, padded, incl. cover, adjustable angle, depth and width	Thigh guides with knee-support for size 1 12 cm long, padded, incl. cover, adjustable angle, depth and width	Thigh guides with knee-support for size 2 17 cm long, padded, incl. cover, adjustable angle, depth and width
Price	Price	Price
☐ 73 00 049 £ 134	☐ 73 01 049 £ 134	☐ 73 02 049 £ 134


Configuration option - seat

Standard seat pad – incl. cover (standard configuration, if nothing else stated) 	Seat pad – incl. cover with anatomic leg guidance 	Adjustable pelvic guidance padded, incl. cover 	
Price	Upcharge		Price
☐ for size 0 (73 00 018) incl. ☐ for size 1 (73 01 018) incl. ☐ for size 1 b (73 01 118) incl. ☐ for size 2 (73 02 018) incl.	☐ for size 0 (73 00 019) £ 77 ☐ for size 1 (73 01 019) £ 77 ☐ for size 1 b (73 01 119) £ 87 ☐ for size 2 (73 02 019) £ 87	☐ for size 0 – seat width A = 15 – 24 cm (73 00 021) incl. ☐ for size 1 + 1 b – seat width A = 18 – 28 cm (73 01 021) incl. ☐ for size 2 – seat width A = 26 – 36 cm (73 02 021) incl.	

Seat option

Leg abductors – including leg length compensation, padded, incl. cover
(already includes thigh guides)

	for size 0	for size 1+1b	for size 2
A	0° – 20°	0° – 20°	0° – 20°
B	7 – 10 cm	9 – 12 cm	10 – 14 cm
C	3 cm	3 cm	3 cm
	Price		
☐ for size 0 (73 00 022)	£ 147		
☐ for size 1 (73 01 022)	£ 147		
☐ for size 1 b (73 01 122)	£ 147		
☐ for size 2 (73 02 022)	£ 147		



Accessories

Pommel – quick release, padded incl. cover

	Width	Depth
Size 0	7 cm	10 cm
Size 1	8 cm	14 cm
Size 2	10 cm	17 cm
	Price	
☐ Size 0 (73 00 023)	£ 99	
☐ Size 1 (73 01 023)	£ 99	
☐ Size 2 (73 02 023)	£ 99	

Split footplates

Split footplate – individually adjustable, adjustable in height, width and
angle, foldable, fixed to the seat unit, with non-slip surface and 4 slots
(each) for attaching straps

	Upcharge / pair
☐ for size 0 (73 00 052)	£ 142
☐ for size 1 (73 01 052)	£ 142
☐ for size 1 b (73 01 152)	£ 142
☐ for size 2 (73 02 052)	£ 142

Split footplate with sandals – individually adjustable in height, width,
angle and flip away

	Upcharge / pair
☐ Size 0 (73 00 051)	£ 163
☐ Size 1 (73 01 051)	£ 163
☐ Size 2 (73 02 051)	£ 163


Accessories

Therapy tray – made of wood, with edge, swing-away 			Therapy tray – made of wood, with edge, fixed on both sides with quick release 			Therapy tray – transparent, with edge, swing-away 		
	Width	Depth		Width	Depth		Width	Depth
for size 0+1	50 cm	40 cm	for size 0+1	50 cm	40 cm	for size 0+1	50 cm	40 cm
for size 2	60 cm	50 cm	for size 2	60 cm	50 cm	for size 2	60 cm	50 cm
Price			Price			Upcharge / set		
☐ for size 0 (73 00 024)	£ 134		☐ for size 0 (73 00 050)	£ 150		☐ for size 0 (73 00 037)	£ 178	
☐ for size 1 (73 01 024)	£ 134		☐ for size 1 (73 01 050)	£ 150		☐ for size 1 (73 01 037)	£ 178	
☐ for size 2 (73 02 024)	£ 134		☐ for size 2 (73 02 050)	£ 150		☐ for size 2 (73 02 037)	£ 178	



Accessories

4 Castors – 2 with locking devices 		4 Castors – all with locking device 		Central parking brake for rear castors incl. 4 dual castors (100 mm), operated by using the parking lever which is mounted to the seat system  	
Upcharge / set		Upcharge / set		Price / set	
☐ 75 mm (73 00 058) incl.		☐ 75 mm (73 00 060)	£ 22	☐ for size 0+1 parking lever mounted left	(74 01 051) £ 210
☐ 100 mm (73 01 058)	£ 37	☐ 100 mm (73 01 060)	£ 58	☐ for size 0+1 parking lever mounted right	(74 01 551) £ 210
☐ 125 mm (73 02 058)	£ 58	☐ 125 mm (73 02 060)	£ 78	☐ for size 1b+2 parking lever mounted left	(74 02 051) £ 210
				☐ for size 1b+2 parking lever mounted right	(74 02 551) £ 210
				☐ for size 0+1 parking lever mounted left – castors front with brakes	(74 01 151) £ 240
				☐ for size 0+1 parking lever mounted right – castors front with brakes	(74 01 651) £ 240
				☐ for size 1b+2 parking lever mounted left – castors front with brakes	(74 02 151) £ 240
				☐ for size 1b+2 parking lever mounted right – castors front with brakes	(74 02 651) £ 240

Accessories

Hydraulic height adjustment (retrofitting not possible) 		Electric height adjustment by a 24-volt electric motor 		Additional gas dampers (2 pcs) to support the lift function, not possible together with hydraulic height adjustment! 2x 300 N recommended from 30 kg load 2x 500 N recommended from 40 kg load 		Anti tip – adjustable in height 	
Upcharge / set		Price / set		Price / set		Price	
☐ for size 0+1 (74 01 050)	£ 204	☐ 74 00 002	£ 458	☐ 2x 300 N (74 00 003)	£ 75	☐ 74 00 004 (for standard & electr. height adjustment)	£ 99
☐ for size 1b+2 (74 02 050)	£ 204			☐ 2x 500 N (74 01 003)	£ 75	☐ 74 01 004 (for hydr. height adjustment)	£ 99

Accessories

Wheel base extension inserts at front, length = 10 cm 		Incontinence coating in the seat area 		Harness assembly kit for the seat consisting of 2 brackets with clip buckle and 2 clip buckles (drilled) 		Harness assembly kit for the back consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled) 	
Price / set		Price		Price / set		Price / set	
☐ 74 00 005	£ 95	☐ 73 00 020	£ 66	☐ 20 mm (73 00 025)	incl.	☐ 20 mm (73 00 026)	incl.
				☐ 25 mm (73 01 025)	incl.	☐ 25 mm (73 01 026)	incl.

Padding

Padded cover for footplate – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps) 		Padded cover for split footplate – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps) 		Padded tray cover – covered with black Flexi cover, height: 1 cm 	
Price		Price		Price	
☐ for size 0 (73 00 038)	£ 68	☐ for size 0 (73 00 039)	£ 89	☐ for size 0 (73 00 054)	£ 78
☐ for size 1 (73 01 038)	£ 68	☐ for size 1 (73 01 039)	£ 89	☐ for size 1 (73 01 054)	£ 78
☐ for size 2 (73 02 038)	£ 68	☐ for size 2 (73 02 039)	£ 89	☐ for size 2 (73 02 054)	£ 78


Harness

The **madita-fun.** harness is assembled with the harness assembly sets (see page 12). If a colour is not selected, we will use black as the standard colour for your order.

Harness

Chest harness – made of elastic, breathable material, three straps with clips, sewn strap to the left shoulder


	Size 0	Size 1	Size 2	Size 3
A	13 cm	13 cm	18 cm	18 cm
B	17.5 cm	19 cm	20 cm	21 cm
C	15 cm	15 cm	19 cm	28 cm
D	3.5 cm	5 cm	5 cm	5 cm
E	25.5 cm	28 cm	34 cm	44 cm
Strap width	2 cm	2 cm	2.5 cm	2.5 cm



	Price	Colour
☐ Size 0 (73 00 032)	£ 104	☐ Black (91 00 040)
☐ Size 1 (73 01 032)	£ 104	☐ Light grey (91 00 049)
☐ Size 2 (73 02 032)	£ 104	
☐ Size 3 (73 03 032)	£ 104	

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips

	Size 0	Size 1	Size 2
A	11 cm	11 cm	18 cm
B	21 cm	21 cm	22 cm
C	12 cm	12 cm	15 cm
D	2.5 cm	3.5 cm	3.5 cm
E	25.5 cm	28 cm	34 cm
Strap width	2 cm	2 cm	2.5 cm



	Price	Colour
☐ Size 0 (73 00 056)	£ 119	☐ Black (91 00 040)
☐ Size 1 (73 01 056)	£ 119	☐ Light grey (91 00 049)
☐ Size 2 (73 02 056)	£ 119	

Harness

T-shape groin harness – made of elastic, breathable material, with two clips at front, straps sewn at back


	Size 0	Size 1	Size 2
A	23 cm	25 cm	28 cm
B	15 cm	17 cm	20 cm
C	33 cm	35 cm	38 cm
D	24 cm	26 cm	28 cm
Strap width	2 cm	2 cm	2.5 cm



	Price	Colour
☐ Size 0 (73 00 033)	£ 89	☐ Black (91 00 040)
☐ Size 1 (73 01 033)	£ 89	☐ Light grey (91 00 049)
☐ Size 2 (73 02 033)	£ 89	


Groin harness (standard) – made of elastic, breathable material, with two clips at front, straps sewn at back

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2.5 cm	2.5 cm



	Price	Colour
☐ Size 1 (73 01 034)	£ 89	☐ Black (91 00 040)
☐ Size 2 (73 02 034)	£ 89	☐ Light grey (91 00 049)
☐ Size 3 (73 03 034)	£ 89	

4-point pelvic harness – with pad and adjustable clip on both sides




	Size 0	Size 1	Size 2
Pad width	4 cm	4 cm	5 cm
Pad length per side	16 cm	18 cm	28 cm
Strap length per side	30 cm	30 cm	40 cm
Strap width	2 cm	2.5 cm	2.5 cm

	Price	Colour
☐ Size 0 (73 00 036/UK)	incl.	☐ Black (91 00 040)
☐ Size 1 (73 01 036/UK)	incl.	☐ Light grey (91 00 049)
☐ Size 2 (73 02 036/UK)	incl.	


Harness

Foot cuff – with deflection pulley and Velcro®, colour: black (not possible in combination with padded cover for footplate/split footplate)



	Ankle circumference	Width Gurt
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm
	Price / pair	
☐ Size 1 (73 01 061)	£ 85	
☐ Size 2 (73 02 061)	£ 85	
☐ Size 3 (73 03 061)	£ 85	

Foot positioner – with straps, incl. velcro fastener (not possible in combination with padded cover for footplate/split footplate)



	A	Strap width
Size 1	8 cm	2 cm
Size 2	10 cm	2 cm
Size 3	13 cm	2.5 cm
	Price / pair	
☐ Size 1 (73 01 062)	£ 75	
☐ Size 2 (73 02 062)	£ 75	
☐ Size 3 (73 03 062)	£ 75	

madita-fun mini.

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____



Signature: _____

Reference:
The therapy chair for initial assistance.

- incl. armrests, footplate with upholstery, ergonomically shaped push handle, standard seat pad incl. cover
- 4 Castors (75 mm) – 2 with locking devices
- incl. back unit low / high in standard design
- incl. pelvis harness
- Height adjustable per gas pressure spring (800N)
- Frame colour: Telegrey 4 (white grey)
- Adjustment possibilities:
 - Seat angle (by gas pressure spring 250N, (-)5° – 37° / 0° – 42°)
 - Back angle
 - Seat depth and seat width
 - Armrest, height, width and angle adjustable
 - Back, height adjustable
 - Footplate, knee angle and height adjustable


Dimensions - basic model


A	Seat depth	12 – 19 cm
B	Seat width	14 – 22 cm
C	Back height with short back	28 – 33 cm
C	Back height with long back	30 – 35 cm
D	Seat height	25 – 62 cm
E	Lower leg length	12 – 22 cm
F	Total width	52 cm
G	Total length	58 cm
	Tilt-in-space	(-)5° – 37° / 0° – 42°
	Back angle	(-)5° – 25°
	Max. load	20 kg
	Weight	16 kg



Basic model	
madita-fun mini.	
	Price
☐ 3100 000/UK	£ 1,420


Seat angle - default setting			
Please select a default setting ((-)5° - 37° = Standard configuration, if nothing else stated)			
	Price		Price
☐ (-)5° - 37° (3100 001)	incl.	☐ 0° - 42° (3100 101)	incl.

Covers			
	Pink (No. 114)		Light blue (No. 106)
	Price		Price
☐ 9100 114	incl.	☐ 9100 106	incl.
			Jeans (No. 116)
			Price
		☐ 9100 116	incl.



Back incl. pad and cover
standard

	Height	Width
Back (short)	28 – 33 cm	24 cm
Back (long)	30 – 35 cm	24 cm
		Price
☐ Back (short)	(3100002)	incl.
☐ Back (long)	(3100003)	incl.



Back with integrated upper arm guidance incl. pad and cover – with adjustable hollow form in the lower trunk area through upholstery wedges (height approx. 3cm)

	Height
Back (short)	28 – 33 cm
Back (long)	30 – 35 cm
Price	
☐ Back (short)	(3100005) £ 32
☐ Back (long)	(3100006) £ 32



Lateral pads – padded and covered, incl. bracket


Height	7 cm
Depth	9 cm
Chest width	14 – 20 cm
Price / pair	
☐ 3100009	£ 69





Lateral pads – long (trunk guidance) padded and covered, incl. bracket

Height	14 cm
Depth	6 cm
Chest width	14 – 20 cm
Price / pair	
☐ 3100007	£ 103

Name embroidery on the back pad



please fill out the form on page 25! Embroidery of flexi fabric is not available.

Configuration possibilities - seat			
Seat pad – standard, incl. cover	Seat pad – with anatomic leg guidance, incl. cover		Abduction seat – with a seat plate (up to 10 cm abductable) in the front sitting area, internal and external leg guidance (approx. 30° abduction per side), incl. pad and cover and a 10 cm wider adjustable footplate
			Upcharge
☐ 3100 011	incl.	☐ 3100 012	£ 79
☐ 3100 013	£ 147		
Thigh guides – padded and covered, angle, width, depth adjustable (included in abduction seat)	Internal guide – abductable (included in abduction seat)	Incontinence coating – in the seat area	
			
	Price / pair	Price / set	Price
☐ 13 cm (3100 016)	£ 103	☐ 3100 015	£ 115
☐ 16 cm (3100 017)	£ 103	☐ 3100 014	£ 69

Accessories

Headrest (standard) padded, height adjustable 	Headrest (standard) padded, height, depth and angle adjustable 	Headrest (shell-shaped) padded, height, depth and angle adjustable 	Therapy tray – made of wood, with edge 	Padded tray cover covered with black Flexi cover, height: 1 cm 
Price	Price	Price	Price	Price
☐ 3100020 £ 131	☐ 3100021 £ 158	☐ 3100124 £ 185	☐ 3100022 £ 138	☐ 3100023 £ 78

Accessories


Harness assembly kit for the seat consisting of 2 brackets with clip buckle and 2 clip buckles (drilled) 	Harness assembly kit for the back consisting of 2 threaded plates with clip buckles and 2 clip buckles (drilled) 	The madita-fun mini. harness is assembled with the harness assembly sets. If a colour is not selected, we will use black as the standard colour for your order.
Price / set	Price / set	
☐ 3100030 incl.	☐ 3100031 incl.	

Fixing aids



Harness

Chest harness – made of elastic, breathable and moisture repellent material, three straps with clips, sewn strap to the left shoulder 			T-shape groin harness – made of elastic, breathable and moisture repellent material, with two clips at front, straps sewn at back 		
	Size mini	Size 0			
A	8.5 cm	13 cm	A	18 cm	23 cm
B	17.5 cm	17.5 cm	B	10 cm	15 cm
C	14 cm	15 cm	C	25 cm	33 cm
D	3.5 cm	3.5 cm	D	19 cm	24 cm
E	24.5 cm	25.5 cm	Strap width	2 cm	2 cm
Strap width	2 cm	2 cm			
	Price	Colour		Price	Colour
☐ Size mini (3105032)	£ 104	☐ Black (9100040)	☐ Size mini (7305033)	£ 89	☐ Black (9100040)
☐ Size 0 (3100032)	£ 104	☐ Light grey (9100049)	☐ Size 0 (7300033)	£ 89	☐ Light grey (9100049)

Harness

Pelvic harness – with pad and adjustable clip on both sides 			4-point pelvic harness – with pad and adjustable clip on both sides 		
	Size mini	Size 0		Size mini	Size 0
Pad width	3.5 cm	4 cm	Pad width	3.5 cm	4 cm
Pad length per side	10 cm	16 cm	Pad length per side	10 cm	16 cm
Strap length per side	20 cm	30 cm	Strap length per side	20 cm	30 cm
Strap width	2 cm	2 cm	Strap width	2 cm	2 cm
	Price	Colour		Upcharge	Colour
☐ Size mini (7305035)	incl.	☐ Black (9100040)	☐ Size mini (7305036)	£ 32	☐ Black (9100040)
☐ Size 0 (7300035/UK)	incl.	☐ Light grey (9100049)	☐ Size 0 (7300036/UK)	£ 32	☐ Light grey (9100049)

Harness

Foot cuff – with deflection pulley and Velcro®, colour: black (not possible in combination with padded cover for footplate) 			Foot fixing – with straps, incl. velcro fastener (not possible in combination with padded cover for footplate) 		
	Ankle circumference	Harness width		A	Strap width
Size 1	12 – 19 cm	25 mm	Size 1	8 cm	2 cm
Size 2	15 – 23 cm	25 mm	Size 2	10 cm	2 cm
Size 3	19 – 27 cm	40 mm	Size 3	13 cm	2.5 cm
	Price / pair			Price / pair	
☐ Size 1 (7301061)	£ 85		☐ Size 1 (7301062)	£ 75	
☐ Size 2 (7302061)	£ 85		☐ Size 2 (7302062)	£ 75	
☐ Size 3 (7303061)	£ 85		☐ Size 3 (7303062)	£ 75	

till.

I am your adjustable supine standing frame.



new

Size
3 & 4

Further information can
be found from page 40.

madita hd.

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____

Reference:
The reliable therapy chair.

- reinforced short back incl. padding and cover
- incl. armrests, footplate and push handle
- incl. pelvic harness
- incl. 25mm harness assembly kit for the seat and back
- incl. 1x weight-bar (ca. 8 kg) and 4 castors (75 mm) with locking devices
- Frame colour: anthracite metallic
- Adjustment options:
 - height adjustment via gas pressure spring (350 N)
 - tilt-in-space
 - back angle
 - seat depth and width
 - angle and width adjustable armrests
 - knee angle, foot angle and height adjustable footplate


Basic model
madita hd.

Price

☐ **30 02 100/UK**

£ 1,590

Covers

**Black/
Dark grey**
(No. 105 + 92)

Price

☐ **91 00 205** incl.

**Dark grey /
Deep orange**
(No. 92 + 107)

Price

☐ **91 00 213** incl.

**Dark grey / Light
blue**
(No. 92 + 106)

Price

☐ **91 00 214** incl.


Coated covers

**Flexi Black/Light
grey**
(No. 43 + 50)

Upcharge

☐ **91 00 226** £ 152

Dimensions - basic model

		madita hd.	
A	Seat depth	30 – 42 cm	
B	Seat width	28 – 38 cm	
C	Back height with short back	42 – 52 cm	
C	Back height with long back	53 – 63 cm	
D	Seat height	46 – 58 cm	
E	Lower leg length	30 – 48 cm	
F	Total width	53 cm	
G	Total length	77 cm	
	Tilt-in-space	(-)5° – 20°	
	Back angle	(-)5° – 25°	
	Max. load	65 kg	
	Weight	29 kg	


Headrest

Headrest

with reinforced bracket, padded incl. cover, height and angle adjustable



Headrest (lateral guide, angle adjustable) – padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable



A = Width contact surface
B = Height Headrest
C = Width lateral guide

new


Price
£ 199

□ 30 02 104

	A	B	C
Size 1	13.5cm	10.5cm	10 cm
Size 2	15cm	10.5cm	13.5cm
Size 3	17 cm	12.5cm	16.5cm
	Price		
□ Size 1 (30 01 030)	£ 348		
□ Size 2 (30 02 030)	£ 348		
□ Size 3 (30 03 030)	£ 348		

Back

Back – reinforced, incl. pad and cover



	Height	Depth
Short back	42 – 52 cm	Size 1 8 cm
Long back	53 – 63 cm	Size 2 11 cm
	Price	Price / pair
□ short (30 02 102) incl.	£ 95	□ Size 1 (30 01 025)
□ long (30 02 103) incl.	£ 95	□ Size 2 (30 02 025)

Lateral pads

Lateral pads – with reinforced holder, padded, incl. cover



Tailored lateral pads
with reinforced bracket,
padded, incl. cover

Height: ____ cm
Depth: ____ cm

Price / pair

□ 30 09 025

£ 168

Lateral pads



Lateral pads size 1, long (trunk guidance)
with reinforced bracket,
padded, incl. cover

Height: 15 cm
Depth: 8 cm

Price / pair

□ 30 01 026

£ 158

Lateral pads size 2, long (trunk guidance)
with reinforced bracket,
padded, incl. cover

Height: 17 cm
Depth: 9 cm

Price / pair

□ 30 02 026

£ 158

Tailored lateral pads, long (trunk guidance)
with reinforced bracket,
padded, incl. cover


Height: ____ cm
Depth: ____ cm

Price / pair

□ 30 09 026

£ 220

Configuration possibilities - seat


Standard seat pad – incl. cover (standard configuration, if nothing else stated)	Seat pad – incl. cover with anatomic leg guidance	Incontinence coating – in the seat area
		
Price	Upcharge	Price
☐ 3002018 incl.	☐ 3002019 £ 85	☐ 3000001 £ 68

Tilt-in-space adjustment

Thigh guides

Trigger lever mounted directly on the gas spring (standard configuration, if nothing else stated – not possible with removable back unit)	Trigger lever mounted on the push bar	Thigh guides padded, incl. cover, adjustable angle, depth and width	
Price	Price	Price / pair	
☐ 7300085 incl.	☐ 7301085 incl.	☐ 17 cm (7302114)	£ 102
		☐ 22 cm (7302115)	£ 102

Accessories

Adjustable pelvic guidance – padded, incl. cover, seat width A = 26 – 36 cm	Pommel – quick release, padded incl. cover	Therapy tray – made of wood with edge, pluggable, fixed on both sides
		
Upcharge	Price	Price
☐ 7302021 £ 209	☐ 7302023 £ 99	☐ 3002023 £ 150

Accessories

Weight ca. 8 kg (max. 3 extra weights possible)	4 Castors – 100 mm, all with locking devices	Wheel base extension length = 10cm, inserts at back
		
Price	Quantity	Upcharge / set
☐ 3002021 £ 77	☐ 7301060	☐ 7301028 £ 98


Padding

Padded footplate cover – covered with black Flexi cover, height: 3 cm (not possible in combination with foot straps)	Padded tray cover – covered with black Flexi cover, height: 1 cm
	
Price	Price
☐ 7302038 £ 68	☐ 7302054 £ 78

Harness

T-shape groin harness – made of elastic, breathable material, with two clips at front, straps sewn at back

	Size 1	Size 2
A	25 cm	28 cm
B	17 cm	20 cm
C	35 cm	38 cm
D	26 cm	28 cm
Strap width	2 cm	2.5 cm



	Price	Colour
☐ Size 1 (73 01 033)	£ 89	☐ Black (91 00 040)
☐ Size 2 (73 02 033)	£ 89	☐ Light grey (91 00 049)

Groin harness (standard) – made of elastic, breathable material, with two clips at front, straps sewn at back

	Size 1	Size 2	Size 3
A	24 cm	29 cm	33 cm
B	14 cm	17 cm	22 cm
C	35 cm	43 cm	49 cm
Strap width	2 cm	2.5 cm	2.5 cm



	Price	Colour
☐ Size 1 (73 01 034)	£ 89	☐ Black (91 00 040)
☐ Size 2 (73 02 034)	£ 89	☐ Light grey (91 00 049)
☐ Size 3 (73 03 034)	£ 89	

Harness


Pelvic harness – with pad and adjustable clip on both sides



	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2.5 cm	2.5 cm

	Price	Colour
☐ Size 1 (73 01 035/UK)	incl.	☐ Black (91 00 040)
☐ Size 2 (73 02 035/UK)	incl.	☐ Light grey (91 00 049)

4-point pelvic harness – with pad and adjustable clip on both sides




	Size 1	Size 2
Pad width	4 cm	5 cm
Pad length per side	18 cm	28 cm
Strap length per side	30 cm	40 cm
Strap width	2.5 cm	2.5 cm

	Upcharge	Colour
☐ Size 1 (73 01 036/UK)	£ 32	☐ Black (91 00 040)
☐ Size 2 (73 02 036/UK)	£ 32	☐ Light grey (91 00 049)

Harness

Chest harness – made of elastic, breathable material, three straps with clips, sewn strap to the left shoulder


	Size 1	Size 2	Size 3
A	13 cm	18 cm	18 cm
B	19 cm	20 cm	21 cm
C	15 cm	19 cm	28 cm
D	5 cm	5 cm	5 cm
E	28 cm	34 cm	44 cm
Strap width	2 cm	2.5 cm	2.5 cm



	Price	Colour
☐ Size 1 (73 01 032)	£ 104	☐ Black (91 00 040)
☐ Size 2 (73 02 032)	£ 104	☐ Light grey (91 00 049)
☐ Size 3 (73 03 032)	£ 104	

Positioning jacket – with zipper, made of elastic, breathable material, shoulder straps sewn at back, pelvis straps with clips


	Size 1	Size 2
A	11 cm	18 cm
B	21 cm	22 cm
C	12 cm	15 cm
D	3.5 cm	3.5 cm
E	28 cm	34 cm
Strap width	2 cm	2.5 cm



	Price	Colour
☐ Size 1 (73 01 056)	£ 119	☐ Black (91 00 040)
☐ Size 2 (73 02 056)	£ 119	☐ Light grey (91 00 049)

Harness


Foot cuff – with deflection pulley and Velcro®, colour: black (not possible in combination with padded cover for footplate)



	Ankle circumference	Width Gurt
Size 1	12 – 19 cm	25 mm
Size 2	15 – 23 cm	25 mm
Size 3	19 – 27 cm	40 mm

	Price / pair
☐ Size 1 (73 01 061)	£ 85
☐ Size 2 (73 02 061)	£ 85
☐ Size 3 (73 03 061)	£ 85

Foot fixing – with straps, incl. velcro fastener (not possible in combination with padded cover for footplate)



	A	Strap width
Size 1	8 cm	2 cm
Size 2	10 cm	2 cm
Size 3	13 cm	2.5 cm

	Price / pair
☐ Size 1 (73 01 062)	£ 75
☐ Size 2 (73 02 062)	£ 75
☐ Size 3 (73 03 062)	£ 75

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____


Signature: _____

Reference:
The cradle subframe with cradle system.

- Electrical height adjustment via hand switch
- 12" inflated rear wheels with drum brake, 8" PU front wheels
- Universal mounting points for all types of cradles
- Frame colour: white aluminium (silver)


Dimensions - cradle subframe

		for three-piece cradle system
A	Height of cradle support from floor	53 - 79 cm
B	Adjustable back angle	0° - 50°
C	Seat angle (adjustable in steps)	0° - 18°
D	Adjustable foot angle	0° - 50°
	Base frame width	64 cm
	Max. load	100 kg
	Weight	44 kg


Basic model cradle subframe


tina. (for three-piece cradle system)
with electrical back adjustment, adjustable seat and foot angle

☐ **4000 101**

Price

£ 2,990

Dimensions - cradle system

		Size 2	Size 3	Size 4	
A	Width	70 cm	70 cm	70 cm	
B	Length of head section	60 cm	75 cm	90 cm	
C	Centre section length	35 cm	35 cm	35 cm	
D	Length of foot section	45 cm	55 cm	65 cm	
E	Total length	140 cm	165 cm	190 cm	
F	Height	12 cm	12 cm	12 cm	

Cradle system

Head section



Note: Size 4 includes an extended push handle and anti-tip device

	Length	Width
Size 2	60 cm	70 cm
Size 3	75 cm	70 cm
Size 4	90 cm	70 cm

Custom modification is possible.

	Price
☐ Size 2 (40 02 102)	£ 259
☐ Size 3 (40 03 102)	£ 259
☐ Size 4 (40 04 102)	£ 548

Seat section



	Length	Width
Size 2-4	35 cm	70 cm

Custom modification is possible.

	Price
☐ Size 2-4 (40 02 103)	£ 179

Leg section



	Length	Width
Size 2	45 cm	70 cm
Size 3	55 cm	70 cm
Size 4	65 cm	70 cm

Custom modification is possible.

	Price
☐ Size 2 (40 02 104)	£ 259
☐ Size 3 (40 03 104)	£ 259
☐ Size 4 (40 04 104)	£ 259

Mattresses

Vacuum mattress – for tina. incl. valve and electric pump with two-way switch



	Length	Width
Size 2	ca. 150 cm	ca. 75 cm
Size 3	ca. 175 cm	ca. 75 cm
Size 4	ca. 200 cm	ca. 80 cm

Cover made of breathable fabric with insert made of spacer fabric (for adjacent vacuum mattress)

Support mattress – with edge of foam



	Length	Width
Size 2	140 cm	70 cm
Size 3	165 cm	70 cm
Size 4	190 cm	70 cm

Custom modification is possible.


	Price		Price		Price
☐ Size 2 (41 02 202)	£ 1,990	☐ for size 2 (41 02 201)	£ 249	☐ Size 2 (40 02 105)	£ 419
☐ Size 3 (41 03 202)	£ 1,990	☐ for size 3 (41 03 201)	£ 249	☐ Size 3 (40 03 105)	£ 419
☐ Size 4 (41 04 202)	£ 1,990	☐ for size 4 (41 04 201)	£ 249	☐ Size 4 (40 04 105)	£ 419

Accessories

Extended push handle – for head sections larger than 75 cm (already included with head section size 4)



Anti tip – pivoting, adjustable in height (already included with head section size 4)



	Upcharge		Price
☐ 40 00 120	£ 99	☐ 40 00 140	£ 190

Accessories

Tray for medical devices – height adjustable, with slotted holes for strap (25 mm) to secure the medical devices



Width	43cm
Depth	40 cm


Infusion holder – telescopic



	Price		Price
☐ 40 09 011	£ 290	☐ 40 00 010	£ 185

Large support cushion


Wide seat belt – strap with sliding unit to fix on the tube frame



Please provide the required total length (A): _____ cm

	Price		Price
☐ 40 00 005	£ 115	☐ 40 00 009	£ 155

Body measurements - seating



The following measurements can be preset:

A Body size	_____	F Knee width	_____
B Sternum height	_____	G Knee height	_____
C Inside leg length	_____	H Seat depth	_____
D Chest width	_____	I Elbow height	_____
E Pelvis width	_____	J Back height	_____

Name embroidery

Our **smilla.**, **madita-fun.** & **madita-fun mini.** seating aids can be customised by embroidering the child's name on the back pad.

Simply choose the name embroidery option and tell us the desired name.

- The height of the embroidery is 2.5 cm - width dependant on name length
- The embroidery is placed on the right side of the back pad
- The embroidery is in a silver-white thread
- The font of the embroidery is shown in the example picture
- **Note: embroidery in not available with Flexi fabric**



Configuration name embroidery

Indication of the name
please in BLOCK LETTERS
(max. 12 letters)

--	--	--	--	--	--	--	--	--	--	--	--	--

carrot 3.

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____

Reference:
The most comfortable and adaptable car seat.


- Anatomically shaped car seat with depth-adjustable headrest, integrated armrests, padded pelvic support, 5-point seat belt with quick-release, safety belt guide left and right
- Incl. headrest wedges with Velcro fastener, seat and backrest upholstery and 2 torso wedges
- Removable and washable covers
- Incl. shoulder module for size L and backrest extension (5 cm) for size XL
- Tested and approved according to the ECE-44/03 standard
- Possible settings :
 - Backrest height +10 cm
 - seat depth +3,5 cm
 - depth-adjustable headrest
 - Sleeping position



Crash-test up to
49 kg for the L
version and 75 kg
for the XL version!

Dimensions - basic model

		Size S	Size L	Size XL
A	Seat depth	25 - 28,5 cm	25 - 28,5 cm	44 cm
B	Seat width	26 cm	26 cm	38 cm
C	Backrest height	30 - 40 cm	45 - 55 cm	50 - 60 cm
D	Headrest height	52 - 62 cm	67 - 77 cm	73,5 - 83,5 cm
E	Total width	42 cm	42 cm	42 cm
	Footrest height	25 - 46 cm		
	Max. load*	10 - 50 kg	10 - 50 kg	40 - 75 kg
	Body height	up to 145 cm	up to 160 cm	up to 176 cm
	Weight**	from 8 kg	from 8 kg	from 8 kg



* The **carrot3.** has been tested to ECE R44 / 04 from 15 to 36 kg. In addition, dummies of 10, 49 and 75 kg were tested with success in accordance with the standard.

** In cases with support from GB or consultant.

Basic model					
carrot 3. – Size S		carrot 3. – Size L incl. shoulder module and headrest L		carrot 3. – Size XL incl. shoulder module, backrest extension (5cm) and headrest L	
		Price		Price	
☐ Red	(BRCAR/1200)	£ 1,249	☐ Red	(BRCAR/1201)	£ 1,449
☐ Black	(BRCAS/1200)	£ 1,249	☐ Black	(BRCAS/1201)	£ 1,449
		Price		Price	
				☐ Black	(BRCAS/1230)
				£ 1,689	

Configuration possibilities - seat					
Rodded Waterfall Seat Base Extension – for sizes S + L		Rodded Seat Extension (5 cm) – for sizes S + L (impossible in combination with the leg rest and foot rest)		Seat Extension – for sizes S + L, insertable, max. extension = 10 cm (can only be used with an insertable seat extension (5 cm), leg rest or foot rest)	
					
		Price		Price	
☐ Red	(BRCAR/1202)	£ 110	☐ Red	(BRCAR/1207)	£ 110
☐ Black	(BRCAS/1202)	£ 110	☐ Black	(BRCAS/1207)	£ 110
		Price		Price	
				☐ Red, 5 cm	(BRCAR/1204)
				☐ Black, 5 cm	(BRCAS/1204)
				☐ Red, 10 cm	(BRCAR/1205)
				☐ Black, 10 cm	(BRCAS/1205)
				£ 72	
				£ 72	
				£ 89	
				£ 89	

Accessoires					
Backrest Extension – for sizes S + L, insertable, max. extension = 10 cm (backrest extension (5 cm) for size XL already included in the basic model)		3D Body Balance Pads –made of foam, for thoracic guidance, pelvic or thigh, with velcro, grey		Air-Mesh Cushions/Pads for little ones	
					
Price		Price / pair		Price	
☐ Red, 5 cm	(BRCAR/1217)	£ 45	☐ BRCASB/0073	£ 59	☐ BRCASB/0067
☐ Black, 5 cm	(BRCAS/1217)	£ 45		£ 79	☐ Red (BRCAR/1206*)
☐ Red, 10 cm	(BRCAR/1218)	£ 59			☐ Black (BRCAS/1206*)
☐ Black, 10 cm	(BRCAS/1218)	£ 59			
		Price		Price	

Accessoires					
Supporting tray – with side guide, to be attached to the car seat with Velcro		Plug-in abduction block – for sizes S + L (INCOMPATIBLE with leg rest, seat extension (5 cm) or foot rest)		Abduction block – insertable, for sizes S + L (COMPATIBLE with leg rest, seat extension (5 cm) or foot rest)	
					
		Price		Price	
☐ Red, for size S+L	(BRCAR/1213)	£ 149	☐ BRCA/1225	£ 195	☐ BRCA/1211
☐ Black, for size S+L	(BRCAS/1213)	£ 149			
☐ Black for size XL	(BRCA/1238)	£ 199			
		Price		Price	

Accessoires

Foot Rest - for sizes S + L, pluggable, swivelling, height-adjustable, with cover and Velcro fastening for legs and feet

	Longueur mollets		
short	25 - 31 cm		
medium	30 - 41 cm		
long	35 - 46 cm		
			Price
☐ short	(BRCAG/1208)		£ 180
☐ medium	(BRCAG/1209)		£ 180
☐ long	(BRCAG/1210)		£ 180

Attachments ISOFIX

Swivel base with ISOFIX attachment - for sizes S + L, with stepless backrest angle adjustment and reinforcement foot (adjustable from 34 to 40 cm), with cushion to protect the car seat



Extended brace for ISO-turn adjustable from 36 to 55 cm

Swivel Base Footplate Adaptor



ISO-Base - for mounting the **carrot3.** in all sizes, for mounting on the driver's and passenger's side, without swivel function



Price / set			Supplement		Price		Price	
☐ Driver's side	(BRCSB/0005/F)	£ 599	☐ BRCBSB/0054	£ 119	☐ BRCBSB/0037	£ 99	☐ BRCBSB/0007	£ 235
☐ Passenger side	(BRCBSB/0005/B)	£ 599						

Growing kits

Growth Kit (from size S to size L) - to convert a **carrot 3.** (size S) to size L, consisting of a shoulder module, headrest (large) and thorax protector (large)



Growth Kit (size L to size XL) - to convert a **carrot 3.** (size L) to size XL, consisting of a wide and deep seat, long connecting tubes for the backrest, straps to be attached to the backrest and a harness for the user.
Please remove the existing orange type plate (ECE-44/03) at the rear during conversion, as the application range of 40 to 75 kg is out of standard!



Price / set			Price / set		
☐ Red	(BRCBSR/0070)	£ 695	☐ Red	(BRCAR/1234)	£ 750
☐ Black	(BRCBSB/0070)	£ 695	☐ Black	(BRCAS/1234)	£ 750

Incontinence covers

Incontinence cover for seat and back Dartex black	Incontinence cover for shoulder protection wing set Dartex black	Incontinence cover for chest shield Dartex black
Price	Price	Price
☐ for size S+L (CAR0066) ☐ for size XL (CAR0080)	☐ BRCSB/0056	☐ for size S (CAR0068) ☐ for size XL (CAR0081)
£ 130 £ 130	£ 48	£ 35 £ 35

Incontinence covers

Incontinence cover for headrest Dartex black	Incontinence cover for seat extention Dartex black	Incontinence cover for back extention Dartex black
Price	Price	Price
☐ for size S (CAR0069) ☐ for size XL (CAR0082)	☐ 50mm (BRCSB/0061) ☐ 100mm (BRCSB/0062)	☐ 50mm (BRCSB/0063) ☐ 100mm (BRCSB/0064)
£ 35 £ 35	£ 35 £ 35	£ 35 £ 35

Anti-Escape-Vests

	Anti-Escape-Vest Size S Chest circumference = 60-70 cm	Anti-Escape-Vest Size M Chest circumference = 70-80 cm	Anti-Escape-Vest Size L Chest circumference = 80-90 cm
	Price	Price	Price
	☐ BRCSB/0136S	☐ BRCSB/0136M	☐ BRCSB/0136L
	£ 110	£ 110	£ 110

Accessoires

Covering for security buckle - incl. emergency belt knife, press closure for chest shield, belt



Price
☐ CAR0079
£ 135

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____


Date: _____

Signature: _____

Reference:
The compact and versatile posterior walker.

- Lightweight and folding aluminium frame (size 0 cannot be folded)
- Individual height and width adjustable handles
- EVA rear wheels with integrated backwards stop function in the hub, can be switched off
- EVA front wheels with direction locks
- Friction brake


Dimensions - basic model

		Size 0 frame cannot be folded	Size 1	Size 2	Size 3	Size 4	
A1	Handle height	29 - 53 cm	34 - 58 cm	51 - 75 cm	68 - 94 cm	87 - 113 cm	
A2	Height of forearm supports	-	46 - 70 cm	63 - 87 cm	79 - 103 cm	98 - 124 cm	
B	Handle spacing	22 - 36 cm	22 - 36 cm	30 - 44 cm	38 - 52 cm	43 - 57 cm	
C	Seat height	-	33 cm	45 cm	60 cm	70 cm	
D	Total width	52 cm	51 cm	61 cm	68 cm	75 cm	
E	Total length	58 cm	58 cm	70 cm	77 cm	87 cm	
	Frame height (folded)	37 cm	25 cm	32 cm	34 cm	37 cm	
	Transport dimensions (L x W x H)	58 x 52 x 37 cm	58 x 51 x 25 cm	70 x 61 x 32 cm	77 x 68 x 34 cm	87 x 75 x 37 cm	
	Max. load	20 kg	30 kg	50 kg	75 kg	100 kg	
	Weight	3,9 kg	7,1 kg	8,3 kg	8,5 kg	10,2 kg	

Basic model									
malte mini. size 0 Frame colour: Lazurite blue metallic (frame cannot be folded up)		malte. size 1		malte. size 2		malte. size 3		malte. size 4	
Price		Price		Price		Price		Price	
☐ 35 00 100/UK	£ 470	☐ 75 01 000/UK	£ 495	☐ 75 02 000/UK	£ 495	☐ 75 03 000/UK	£ 515	☐ 75 04 000/UK	£ 570

Frame colours									
Colour selection for size 1			Colour selection for size 2			Colour selection for size 3+4			Special colour (on request) Please provide the requested colour here: RAL _____
									
Price			Price			Price			Price
☐ Telemagenta (pink) (92 00 012)	incl.		☐ Lazurite blue metallic (92 00 010)	incl.		☐ Anthracite metallic (92 00 014)	incl.		☐ 75 00 090
☐ Sky blue (92 00 013)	incl.		☐ Traffic red (92 00 007)	incl.		☐ Pure white (92 00 100)	incl.		£ 90

Accessories									
Forearm armrests with handles – adjustable in depth, height and angle – handles adjustable in depth and angle (for size 0 on request)			Forearm armrests with ball grips – adjustable in depth, height and angle – handles adjustable in depth and angle (for size 0 on request)			Grip bar – pivoting, as a cross-bar between the handles (for sizes 1,2 & 3, not available in conjunction with forearm armrests)			
									
	Length	Width		Length	Width				
Size 1	18 cm	6.5 cm	Size 1	18 cm	6.5 cm				
Size 2	25 cm	8 cm	Size 2	25 cm	8 cm				
Price / set			Price / set			Price			
☐ Size 1 (75 01 001)	£ 160		☐ Size 1 (75 01 002)	£ 160		☐ for size 1 (75 01 003)	£ 155		
☐ Size 2 (75 02 001)	£ 160		☐ Size 2 (75 02 002)	£ 160		☐ for size 2 (75 02 003)	£ 155		
						☐ for size 3 (75 03 003)	£ 155		

Accessories					
Handle stops		Combined pelvic and buttock pad for sizes 1-4 – height, width and depth adjustable, PU padded		Combined pelvic and buttock pad for size 0 – height, width and depth adjustable, PU padded	
					
		Width		Width	
		Size 1	20 – 28 cm	for size 0	20 – 28 cm
		Size 2	25 – 33 cm		
		Size 3	30 – 36 cm		
Price / set		Price / set		Price / set	
☐ 75 00 004	£ 52	☐ Size 1 (75 01 005)	£ 199	☐ for size 0 (35 01 103)	£ 199
		☐ Size 2 (75 02 005)	£ 199		
		☐ Size 3 (75 03 005)	£ 199		

Accessories

Handbrake (not available for size 0)  Assembly in combination with a forearm rest (not available in combination with ball grips) 	Handbrake – with extra small brake handles at reduced grip width with separate locking function 	Friction brake 	Push handle  new
Price / set	Price / set	Price / set	Price
☐ for size 1 (75 01008) £ 128 ☐ for size 2 (75 02008) £ 128 ☐ for size 3 (75 03008) £ 128 ☐ for size 4 (75 04008) £ 128	☐ for size 1 (75 01015) £ 255 ☐ for size 2 (75 02015) £ 255 ☐ for size 3 (75 03015) £ 255 ☐ for size 4 (75 04015) £ 255	☐ for size 0+1 (75 00009) incl. ☐ for size 2-4 (75 01009) incl.	☐ for right (75 00017) £ 118 ☐ for left (75 01017) £ 118

Accessories

Anti tip – incl. support roller (for size 0 without support roller) 	Deflector set (not possible for size 0) 	Bumper wheels – for size 0 	Storage basket 
Price / set	Price / set	Price / set	Price
☐ for size 0 (35 00005) £ 62 ☐ for size 1 (75 01010) £ 62 ☐ for size 2-4 (75 02010) £ 62	☐ 75 00012 £ 99	☐ 35 00010 £ 99	☐ for size 0 (35 00011) £ 42 ☐ for size 1-4 (75 00011) £ 42

Fabric collection.

for our seat unit madita-fun.



Dark grey/Pink-Telemagenta
(No. 92 + 114)



Dark grey/Azur
(No. 92 + 110)



Dark grey/Deep orange
(No. 92 + 107)



Dark grey/Light blue
(No. 92 + 106)



Black/Dark grey
(No. 105 + 92)



Jeans
(No. 116)



Flexi Black/Light grey
(No. 43 + 50)



Flexi Black/Red
(No. 43 + 45)



Flexi Blue
(No. 44)

The fabric combinations shown for our covers can be selected directly from our aid **madita-fun**. (page 14).

malte outdoor.

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____


Signature: _____

Reference:
The compact and versatile posterior-walker for outdoor use.

- Lightweight and folding aluminium frame
- Individual height and width adjustable handles
- Drum brake incl. parking handbrake
- EVA back wheels
- EVA front wheels with direction locks


Dimensions - basic model

		Size 1	Size 2	Size 3
A	Handle height	34 - 58 cm	51 - 75 cm	68 - 94 cm
B	Handle spacing	22 - 36 cm	30 - 44 cm	38 - 52 cm
C	Seat height	33 cm	45 cm	60 cm
D	Total width	56 cm	66 cm	73 cm
E	Total length	66 cm	75 cm	82 cm
	Frame height (folded)	25 cm	32 cm	34 cm
	Transport dimensions (LxWxH)	66x56x25 cm	75x66x32 cm	82x73x34 cm
	Max. load	30 kg	50 kg	70 kg
	Weight	7.1 kg	8.3 kg	9.8 kg



Basic model

malte outdoor. size 1 6" tyres (front) and 8" (rear)	malte outdoor. size 2 8" tyres (front) and 10" (rear)	malte outdoor. size 3 8" tyres (front) and 10" (rear)
	Price	Price
☐ 75 01 100	£ 550	☐ 75 02 100 £ 550
		☐ 75 03 100 £ 595

Frame colours

Colour selection for size 1	Colour selection for size 2	Colour selection for size 3	Special colour (on request) Please provide the requested colour here: RAL _____
			
Price	Price	Price	Price
☐ Telemagenta (Pink) (92 00 012) incl. ☐ Sky blue (92 00 013) incl.	☐ Lazurite blue metallic (92 00 010) incl. ☐ Traffic red (92 00 007) incl.	☐ Anthracite metallic (92 00 014) incl. ☐ Pure white (92 00 100) incl.	☐ 75 00 090 £ 90

Accessories

Handle stops	Combined pelvic and buttock support pad – height, width and depth adjustable, PU padded	Folding seat – padded and covered	Handbrake – with extra small brake handles at reduced grip width with separate locking function
			
	Width		
	Size 1	20 - 28 cm	
	Size 2	25 - 33 cm	
	Size 3	30 - 36 cm	
Price / set	Price / set	Price	Upcharge / set
☐ 75 00 004 £ 52	☐ Size 1 (75 01 005) £ 199 ☐ Size 2 (75 02 005) £ 199 ☐ Size 3 (75 03 005) £ 199	☐ for size 1 (75 01 007) £ 144 ☐ for size 2 (75 02 007) £ 144 ☐ for size 3 (75 03 007) £ 144	☐ for size 1 (75 01 016) £ 99 ☐ for size 2 (75 02 016) £ 99 ☐ for size 3 (75 03 016) £ 99

Accessories

Forearm armrests with handles adjustable in depth, height and angle, handles adjusta- ble in depth and angle	Anti tip – incl. support roller	Backwards stop (retrofitting not possible)	Push handle	Storage basket
				
Length Width				
Size 1 18 cm 6,5 cm				
Size 2 25 cm 8 cm				
Price / set	Price / set	Upcharge / pair	Price	Price
☐ Size 1 (75 01 001) £ 155 ☐ Size 2 (75 02 001) £ 155	☐ 75 00 010 £ 62	☐ for size 1 (75 01 013) £ 128 ☐ for size 2+3 (75 02 013) £ 128	☐ for right (75 00 017) £ 118 ☐ for left (75 01 017) £ 118	☐ 75 00 011 £ 40

☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____


Customer:
Notes:
Reference:

The anterior-walker for more independence and mobility.

- Lightweight aluminium frame (sizes 2+3 can be folded)
- Individual height and width adjustable handles
- EVA tyres (150 mm for sizes 0+1, 200 mm for sizes 2+3)


Dimensions - basic model

		Size 0	Size 1	Size 2	Size 3
A1	Handle height	30 – 47 cm	38 – 59 cm	53 – 75 cm	63 – 88 cm
A2	Height of forearm supports	41 – 58 cm	49 – 70 cm	63 – 85 cm	73 – 98 cm
B	Handle spacing	18 – 32 cm	22 – 36 cm	30 – 44 cm	36 – 50 cm
C	Seat height	25 cm	30 cm	45 cm	45 cm
D	Total width	44 cm	48 cm	62 cm	71 cm
E	Total length	36.5 cm	42.5 cm	57.5 cm	61 cm
E	Total length (telescopic wheelbase)	39 – 47 cm	45 – 55 cm	60 – 74 cm	65 – 79 cm
	Max. load	20 kg	30 kg	50 kg	65 kg
	Weight	4.2 kg	5 kg	7.4 kg	9.1 kg



Basic model

marcy. size 0 - 3

	Price
☐ Size 0 (7600000)	£ 380
☐ Size 1 (7601000)	£ 380
☐ Size 2 (7602000)	£ 380
☐ Size 3 (7603000)	£ 380



marcy. size 0 - 3 (telescopic)

Basic features with additional,
adjustable wheelbase

	Price
☐ Size 0 (7600100)	£ 480
☐ Size 1 (7601100)	£ 480
☐ Size 2 (7602100)	£ 480
☐ Size 3 (7603100)	£ 480



Frame colours

Colour selection for sizes 0+1	Colour selection for size 2	Colour selection for size 3	Special colour (on request) Please provide the requested colour here: RAL _____
			
Price	Price	Price	Price
☐ Telemagenta (Pink) (9200012) incl.	☐ Lazurite blue metallic (9200010) incl.	☐ Anthracite metallic (9200014) incl.	☐ 7600090 £ 90
☐ Sky blue (9200013) incl.	☐ Traffic red (9200007) incl.	☐ Pure white (9200100) incl.	

Accessories

Forearm armrests with handles adjustable in depth, height and angle - handles adjustable in depth and angle			Forearm armrests with grip balls adjustable in depth, height and angle - handles adjustable in depth and angle			Seat – fixed, padded and covered incl. back strap			Handbrake		
											
	Length	Width		Length	Width						
Size 1	18 cm	6,5 cm	Size 1	18 cm	6,5 cm						
Size 2	25 cm	8 cm	Size 2	25 cm	8 cm						
Price / set			Price / set			Price			Price / set		
☐ Size 1 (7601001)	£ 160		☐ Size 1 (7601002)	£ 160		☐ for size 0 (7600004)	£ 95		☐ for size 0 (7600005)	£ 125	
☐ Size 2 (7602001)	£ 160		☐ Size 2 (7602002)	£ 160		☐ for size 1 (7601004)	£ 95		☐ for size 1 (7601005)	£ 125	
						☐ for size 2 (7602004)	£ 95		☐ for size 2 (7602005)	£ 125	
						☐ for size 3 (7603004)	£ 95		☐ for size 3 (7603005)	£ 125	

Assembly in combination with the forearm armrest (not possible in combination with grip balls)




Accessories

Handbrake – with extra small brake handles at reduced grip width with separate locking function			Friction brake			Backwards stop		
								
Upcharge / set			Price / set			Upcharge / pair		
☐ for size 0 (7600010)	£ 245		☐ for size 0+1 (7600006)	£ 35		☐ for size 0+1 (7600009)	£ 130	
☐ for size 1 (7601010)	£ 245		☐ for size 2+3 (7602006)	£ 35		☐ for size 2+3 (7602009)	£ 130	
☐ for size 2 (7602010)	£ 245							
☐ for size 3 (7603010)	£ 245							

Accessories

Direction locks			Storage basket			Push handle		
								
Price / pair			Price			Price		
☐ for size 0+1 (7600008)	£ 50		☐ for size 0 (7600007)	£ 39		☐ for right (7600011)	£ 118	
☐ for size 2+3 (7602008)	£ 50		☐ for size 1 (7601007)	£ 39		☐ for left (7601011)	£ 118	
			☐ for size 2+3 (7602007)	£ 39				

Body measurements - walking



The following measurements can be preset:

A Body size	_____	F Knee width	_____
B Sternum height	_____	G Knee height	_____
C Inside leg length	_____	H Seat depth	_____
D Chest width	_____	I Elbow height	_____
E Pelvis width	_____	J Back height	_____



todd.

The moveable standing trainer.

Further information can be found from page 46.



Customer:

Notes:

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____


Signature: _____

Reference:


The height adjustable supine standing frame.

- Height and angle adjustment (0–90°) via gas pressure spring for sizes 1+2
- Height and angle adjustment (0–90°) via electric power for sizes 3+4
(incl. safety shut-off and manual operation)
- Frame shaped to work with mobile hoists
- Upper body support adjustable in length, incl. padding elements
to adjust the length of the upper body support
- Incl. chest strap and pelvic harness
- Adapter for leg rests (standard)
- Leg rests (standard), adjustable in length, with up to
5 cm pelvic rotation adjustment
- Split footrests (standard), which can be lowered to the floor
via height adjustment of the center column
- Base frame with 5 (for size 1–3) or 6 (for size 4) castors (75 mm)
with total locking devices, one of them with directional lock
- Frame colour: Telegrey 4 (white grey), attachment parts: black
- Incl. tool for all settings (adjustments can be
made without tools)
- incl. push handles (size 3+4)

new
Sizes
3 & 4




Dimensions - basic model



		Size 1	Size 2	Size 3	Size 4
A	Body size	75 - 115 cm	100 - 135 cm	120 - 160 cm	145 - 185 cm
B	Upper body support	Hip/thorax width	14 - 27 cm	17 - 31 cm	20 - 37 cm
C		Length	30.5 - 43 cm	38 - 48 cm	49 - 61 cm
D	Leg rests	Foot to knee	20.5 - 30.5 cm	27.5 - 37 cm	32.5 - 45 cm
E		Inside leg to knee	12 - 18.5 cm	16.5 - 24 cm	20 - 29.5 cm
		Knee angle	0° - 65°		
		Hip angle	0° - 45°		
		Abduction (each leg)	0° - 30°		
F		Pelvic rotation adjustment	0 - 5 cm		
G	Feet	Dorsiflexion	0° - 20°		
G		Plantarflexion	0° - 15°		
H		Inversion/Eversion	0° - 15° each direction		
I	Height of lying surface (for transfer)	53.5 - 73.5 cm	60 - 85 cm	71.5 - 93.5 cm	82.5 - 104.5 cm
J	Tilt	0° - 90°			
	Transport dimensions (L x B x H)	95x51x75 cm	105x61x85 cm	125x76x97 cm	152x77x110 cm
	Max. load	30 kg	60 kg	80 kg	100 kg
	Weight	30 kg	35 kg	80 kg	90 kg

Settings



The following measurements should be pre-set on the standing frame:

A Body size	_____	F Pelvic width	_____
B Axilla height	_____	G Hip height	_____
C Shoulder height	_____	H Knee height	_____
D Inside leg height	_____	I Neck width	_____
E Chest width	_____	J Tray height	_____

□ 42 00 020

Price

incl.

Basic model			
till. size 1 for users with a body size of 75 – 115 cm		till. size 2 for users with a body size of 100 – 135 cm	
with gas spring supported height and angle adjustment		with gas spring supported height and angle adjustment	
Price		Price	
☐ 42 01 000	£ 1,650	☐ 42 02 000	£ 1,790

Basic model			
till. size 3 for users with a body size of 120 – 160 cm		till. size 4 for users with a body size of 145 – 185 cm	
with power supported height and angle adjustment		with power supported height and angle adjustment	
with power supported height and angle adjustment and contourable upper body support (divided into shoulder, trunk and pelvis areas – areas can be individually fixed in the desired position with locking levers)		with power supported height and angle adjustment and contourable upper body support (divided into shoulder, trunk and pelvis areas – areas can be individually fixed in the desired position with locking levers)	
Price		Price	
☐ 42 03 000	£ 2,890	☐ 42 03 100	£ 3,390
Price		Price	
☐ 42 04 000	£ 2,990	☐ 42 04 100	£ 3,490


Covers for size 1+2		Covers for size 3+4		Fabric choice
Black/Azur (No. 105 + 110) 	Black/Pink 2 (No. 105 + 115) 	Black/Petrol (No. 105 + 112) 	Black/Dark grey (No. 105 + 92) 	Fabric choice – please enter here the desired fabric out of our fabric collection: Art. No. (see catalogue)
Price	Price	Price	Price	Price
☐ 9100 233	☐ 9100 234	☐ 9100 237	☐ 9100 236	☐ 9100 999
incl.	incl.	incl.	incl.	incl.

Headrests			
Headrest (standard) padded, incl. cover, height, depth and angle adjustment 		Headrest (shell-shaped) padded, incl. cover, height, depth and angle adjustment 	
Price		Price	
☐ for size 1 (42 01 001)	£ 160	☐ for size 1 (42 01 003)	£ 220
☐ for size 2 (42 02 001)	£ 160	☐ for size 2 (42 02 003)	£ 220
☐ for size 3+4 (42 03 001)	£ 160	☐ for size 3+4 (42 03 003)	£ 220
Price		Price	
☐ for size 1 (42 01 002)	£ 220	☐ for size 1 (42 01 002)	£ 220
☐ for size 2 (42 02 002)	£ 220	☐ for size 2 (42 02 002)	£ 220
☐ for size 3+4 (42 03 002)	£ 220	☐ for size 3+4 (42 03 002)	£ 220

Headrests

Headrest (lateral guide, angle adjustable) – padded, incl. cover, height, depth and angle adjustable, side guides are each angle adjustable

Universal adapter for tailored headrests
15x15 mm, height and depth adjustable



A = Width contact surface
B = Height Headrest
C = Width lateral guide


	A	B	C
Size 1	13.5 cm	10.5 cm	10 cm
Size 2	15 cm	10.5 cm	13.5 cm
Size 3	17 cm	12.5 cm	16.5 cm
Price			
☐ Size 1 (42 01 033)	£ 348		
☐ Size 2 (42 02 033)	£ 348		
☐ Size 3 (42 03 033)	£ 348		




☐ **42 01 004** £ 120

Support pads

Lateral pads for lying surface (standard)
incl. upper arm support, flip away, stepless adjustment for height and width, adjustable for rotation



Pelvic pads for lying surface (standard)
incl. forearm rest, flip away, stepless adjustment for height and width, adjustable for rotation



	Height (A)	Depth (C)	suitable for	Chest width		Height (A)	Depth (C)	suitable for	Chest width (B)
Size 1	10.5 cm	9 cm	till. size 1	14 – 27 cm	Size 1	12 cm	9 cm	till. size 1	14 – 27 cm
			till. size 2	17 – 31 cm				till. size 2	17 – 31 cm
Size 2	11.5 cm	11 cm	till. size 2	17 – 31 cm	Size 2	14 cm	10.5 cm	till. size 2	17 – 31 cm
			till. size 3	20 – 37 cm				till. size 3	20 – 37 cm
Size 3	12.5 cm	16 cm	till. size 3	20 – 37 cm	Size 3	17 cm	15.5 cm	till. size 3	20 – 37 cm
			till. size 4	29 – 46 cm				till. size 4	29 – 46 cm
Size 4	13.5 cm	18 cm	till. size 3	20 – 37 cm	Size 4	21 cm	17.5 cm	till. size 3	20 – 37 cm
			till. size 4	27 – 46 cm				till. size 4	27 – 46 cm
Price / pair					Price / pair				
☐ Size 1 (42 01 005)	£ 195				☐ Size 1 (42 01 007)	£ 195			
☐ Size 2 (42 02 005)	£ 195				☐ Size 2 (42 02 007)	£ 195			
☐ Size 3 (42 03 005)	£ 195				☐ Size 3 (42 03 007)	£ 195			
☐ Size 4 (42 04 005)	£ 195				☐ Size 4 (42 04 007)	£ 195			

Support pads



Knee support pads – individual three-dimensional adjustment, with patella recess and positioning aid



Knee support pads – individual three-dimensional adjustment, with adjustable padded support brackets and positioning aid



Padding – for additional support of the knee, to prevent hyperextension

Inside width (A)			Inside width (A)			Length		Ø	
Size 1	8 cm		Size 1	8 cm		Size 1	8 cm		
Size 2	10 cm		Size 2	10 cm		Size 2	8 cm		
Size 3	12 cm		Size 3	12 cm		Size 3	10 cm		
Size 4			Size 3	12 cm		Size 4	12 cm		
Price / pair			Upcharge / pair			Price			
☐ Size 1 (42 01 009)	£ 320	☐ Size 1 (42 01 010)	£ 360	☐ Size 1 (42 01 011)	£ 98				
☐ Size 2 (42 02 009)	£ 320	☐ Size 2 (42 02 010)	£ 360	☐ Size 2 (42 02 011)	£ 98				
☐ Size 3 (42 03 009)	£ 320	☐ Size 3 (42 03 010)	£ 360	☐ Size 3 (42 03 011)	£ 98				
				☐ Size 4 (42 04 011)	£ 98				

Accessories (Note: These articles cannot be retrofitted!)

for size 1+2



for size 3+4



Push bar

Safety bolt – for release lever, prevents unintentional release of height and angle adjustments (Note: only for size 1+2, not retrofittable!)



	Price		Price / set
☐ for size 1+2 (42 01 012)	£ 190	☐ 42 00 013	£ 45
☐ for size 3+4 (42 03 012)	£ 190		

Leg configuration

Adapter for leg rests (standard)

(standard configuration if nothing else stated)



Adapter for leg rests (abduction)

enables individual and continuous abduction of each leg from 0° - 30° (Note: not retrofittable!)



Leg rests (standard)

length adjustable, with 5 cm pelvic rotation adjustment (standard configuration if nothing else stated)



Leg rests (flex) – adjustable in length and angle, with 5 cm pelvic rotation adjustment (Note: not retrofittable!)



Price		Upcharge		Price / pair		Upcharge / pair	
☐ for size 1+2 (42 01 014)	incl.	☐ for size 1+2 (42 01 015)	£ 495	☐ for size 1 (42 01 016)	incl.	☐ for size 1 (42 01 017)	£ 330
☐ for size 3+4 (42 03 014)	incl.	☐ for size 3+4 (42 03 015)	£ 495	☐ for size 2 (42 02 016)	incl.	☐ for size 2 (42 02 017)	£ 330
				☐ for size 3 (42 03 016)	incl.	☐ for size 3 (42 03 017)	£ 380
				☐ for size 4 (42 04 016)	incl.	☐ for size 4 (42 04 017)	£ 380

Foot configuration

Split footplates (standard)

(standard configuration if nothing else stated | Note: not retrofittable!)



Split footplates (flex) – individually adjustable, dorsiflex (20°), plantarflex (15°) and inversion/eversion (+/-15°) (Note: not retrofittable!)



Price / pair		Upcharge / pair	
☐ for size 1 (42 01 018)	incl.	☐ for size 1 (42 01 019)	£ 150
☐ for size 2 (42 02 018)	incl.	☐ for size 2 (42 02 019)	£ 150
☐ for size 3 (42 03 018)	incl.	☐ for size 3 (42 03 019)	£ 150
☐ for size 4 (42 04 018)	incl.	☐ for size 4 (42 04 019)	£ 150

Accessories

Heel edges – – incl. belt + magnetic lock

(Note: cannot be combined with simple foot strap!)



Simple foot strap

Simple foot strap – with magnetic lock

(Note: cannot be combined with heel edges!)



Price / pair		Price / pair	
☐ for size 1 (42 01 021)	£ 99	☐ for size 1 (42 01 022)	£ 69
☐ for size 2 (42 02 021)	£ 99	☐ for size 2 (42 02 022)	£ 69
☐ for size 3 (42 03 021)	£ 99	☐ for size 3 (42 03 022)	£ 69
☐ for size 4 (42 04 021)	£ 99	☐ for size 4 (42 04 022)	£ 69

Therapy tray



Wood therapy tray with edge
removable, depth, height and angle adjustable

Price	
☐ for size 1 (42 01 024)	£ 195
☐ for size 2 (42 02 024)	£ 195
☐ for size 3 (42 03 024)	£ 195
☐ for size 4 (42 04 024)	£ 195



Padded tray cover
covered with black wipe clean material, height: 1 cm

Price	
☐ for size 1 (42 01 025)	£ 95
☐ for size 2 (42 02 025)	£ 95
☐ for size 3 (42 03 025)	£ 95
☐ for size 4 (42 04 025)	£ 95

Harness

Chest strap – with pad and adjustable clip on both sides (standard configuration if nothing else stated), colour: black



	Size 1	Size 2	Size 3
Pad width	4 cm	6 cm	6 cm
Pad length per side	18 cm	25 cm	30 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2.5 cm	2.5 cm	2.5 cm

Price

☐ Size 1 (42 01 029)	incl.
☐ Size 2 (42 02 029)	incl.
☐ Size 3 (42 03 029)	incl.

Wide Chest strap – padded, colour: black



	Size 1	Size 2	Size 3
Pad width	8 cm	10 cm	12 cm
Pad length per side	10 cm	25 cm	35 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2.5 cm	2.5 cm	2.5 cm


Upcharge / piece

☐ Size 1 (42 01 030)	£ 29
☐ Size 2 (42 02 030)	£ 29
☐ Size 3 (42 03 030)	£ 29

Harness

Positioning jacket – made of elastic, breathable and moisture repellent material, shoulder straps and pelvis straps (25 mm) with clips, colour: black

	A	B	C	D
Size 1	14 cm	26.5 cm	14 cm	33 cm
Size 2	15 cm	29.5 cm	16 cm	36 cm
Size 3	17 cm	35.5 cm	19 cm	41 cm
Size 4	18 cm	38.5 cm	21 cm	43.5 cm



Upcharge / piece

☐ Size 1 (42 01 027)	£ 160
☐ Size 2 (42 02 027)	£ 160
☐ Size 3 (42 03 027)	£ 160
☐ Size 4 (42 04 027)	£ 160

Extension for positioning jacket – for the belt attachment of the positioning jacket

recommended from a neck width of

for size 1	11 cm
for size 2	17 cm
for size 3	21 cm
for size 4	30 cm



Price / set

☐ for size 1+2 42 01 028	£ 39
☐ for size 3+4 42 03 028	£ 39

Harness

Pelvic harness – with pad and adjustable clip on both sides (standard configuration if nothing else stated), colour: black



	Size 1	Size 2	Size 3
for Pelvic width	16 – 25 cm	25 – 40 cm	38 – 46 cm
Pad width	4 cm	6 cm	6 cm
Pad length per side	18 cm	25 cm	30 cm
Strap length per side	25 cm	45 cm	65 cm
Strap width	2,5 cm	2,5 cm	2,5 cm

Price

☐ Size 1 (42 01 031)	incl.
☐ Size 2 (42 02 031)	incl.
☐ Size 3 (42 03 031)	incl.

Pelvic harness – by "Pelvi.Loc", colour: black



	Size 1 (XS)	Size 2 (S)	Size 3 (M)	Size 4 (L)
for Pelvic width	16 – 25 cm	25 – 34 cm	32 – 40 cm	38 – 46 cm
Pad width	6,5 cm	7 cm	7,5 cm	8,5 cm
Pad length per side	16 cm	19 cm	22 cm	24 cm
Strap length per side	25 cm	25 cm	45 cm	45 cm
Strap width	2,5 cm	2,5 cm	2,5 cm	2,5 cm

Upcharge / piece

☐ Size 1 (XS) (42 01 032)	£ 190
☐ Size 2 (S) (42 02 032)	£ 190
☐ Size 3 (M) (42 03 032)	£ 190
☐ Size 4 (L) (42 04 032)	£ 190

Accessories



Rechargeable battery (Li-Ion 2.25 Ah) – for basic models with power supported height adjustment (Note: can only be retrofitted by trained personnel)

	Price
☐ 42 01 026	£ 290

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____

Customer:

Notes:

Reference:


The moveable standing trainer.

- **Basic model (standard)** with height and depth adjustable footplate and four height adjustable support castors (75 mm), rear with springs and lock
- **Basic model (with abduction)** with individual adjustable split footplate, two height adjustable support castors (75 mm) in the front and one support castor (30mm) in the back
- incl. chest support pad, rear support pad and three-dimensionally adjustable knee support pads with patella recess
- Mechanical angle adjustment of the central column from 0 - 15°
- Incl. drum brakes, easy operation via locking lever
- Adjustable camber



Please send us the following measurements for optimum presetting of the standing frame

Settings



The following measurements should be pre-set on the standing frame:


A Body size	_____	H Knee height	_____
B Sternum height	_____	E Chest width	_____
C Elbow height	_____	F Pelvic width	_____
D Crotch height	_____	G Knee width	_____

☐ **49 00 020**

Price

incl.

Dimensions - basic model



		Size 1		Size 2			Size 3		Size 4		
	Body size	75 – 110 cm		100 – 130 cm			120 – 155 cm		150 – 175 cm		
	Wheel size	24"	28"	30"	32"	36"	32"	36"	36"	40"	42"
A	Grip height: footplate – hand rims at top	40 – 48 cm	50 – 58 cm	57 – 65 cm	64 – 69 cm	72 – 80 cm	60 – 68 cm	71 – 79 cm	82 – 90 cm	87 – 95 cm	87 – 95 cm
	Recommended from elbow height	50 cm	60 cm	70 cm	80 cm	92 cm	82 cm	92 cm	92 cm	97 cm	99 cm
B	Width of hand rim at top	0°	–	–	–	–	66 cm	66 cm	69 cm	69 cm	69 cm
		3°	–	–	56 cm	55 cm	62 cm	61 cm	65 cm	64 cm	63 cm
		6°	–	–	53 cm	51 cm	58 cm	56 cm	61 cm	59 cm	58 cm
		9°	45 cm	44 cm	50 cm	47 cm	55 cm	52 cm	57 cm	55 cm	53 cm
		12°	42 cm	41 cm	47 cm	43 cm	–	–	–	–	–
		15°	40 cm	–	–	–	–	–	–	–	–
C	Width of hand rims at bottom (min. – max.)	68 – 74 cm	70 – 73 cm	71 – 83 cm	71 – 83 cm	71 – 84 cm	73 – 83 cm	73 – 85 cm	76 – 86 cm	76 – 86 cm	76 – 90 cm
D	Footplate height – centre of knee support pads	19 – 30 cm		28 – 36 cm			34 – 44 cm		42 – 51 cm		
E	Footplate height – rear support pad support	26 – 42 cm		43 – 58 cm			53 – 69 cm		65 – 80 cm		
	Recommended for inside leg length	30 – 47 cm		48 – 63 cm			58 – 74 cm		70 – 84 cm		
F	Abduction split footplate	0° – 6°		0° – 6°			0° – 30°		0° – 30°		
F	Abduction split footplate (with abducted basic model)	0° – 40°		0° – 40°			–		–		
G	Pelvis dept*	11 – 24 cm		11 – 24 cm			15 – 28,5 cm		15 – 28,5 cm		
	Pelvis dept** (adjustment range of pelvis position via crank max. 6 cm)	11 – 20 cm (+6 cm)		11 – 20 cm (+6 cm)			15 – 27 cm (+6 cm)		15 – 27 cm (+6 cm)		
H	Sternum height (height from footplate to centre of chest support pad)	54 – 80 cm		68 – 101 cm			92 – 123 cm		111 – 140 cm		
I	Tray height from footplate	43 – 66 cm		61 – 79 cm			72 – 96 cm		92 – 110 cm		
	Max. load	30 kg		45 kg			65 kg		80 kg		

* with bracket for rear support pad (standard) ** with bracket for rear support pad (Flip away)

Basic model

todd. size 1 for users with a body size 75 – 110 cm		todd. size 2 for users with a body size 100 – 130 cm		todd. size 3 for users with a body size 120 – 155 cm		todd. size 4 for users with a body size 150 – 175 cm	
Price		Price		Price		Price	
☐ standard	(49 01 000) £ 1,570	☐ standard	(49 02 000) £ 1,650	☐ standard	(49 03 000) £ 1,890	☐ standard	(49 04 000) £ 2,090
☐ with abduction	(49 01 100) £ 1,770	☐ with abduction	(49 02 100) £ 1,850				

Frame colours

Dormant red	Melon yellow	Turquoise blue	Lasurit blue metallic	Anthracite metallic	Pure white
					
Price	Price	 Price	Price	Price	Price
☐ 92 00 003 incl.	☐ 92 00 006 incl.	☐ 92 00 016 incl.	☐ 92 00 010 incl.	☐ 92 00 014 incl.	☐ 92 00 100 incl.

Coated covers		Fabric choice	
Flexi Black		Flexi Light grey	
Price		Price	
☐ 9100043	incl.	☐ 9100050	incl.
		Fabric choice – please enter here the desired fabric out of our fabric collection: Art.-No. (see catalogue) _____	☐ 9100999
		incl.	

Support pads				Accessory				
Chest support pad height adjustable			Chest support pad with lateral guide can be adjusted in height, width and depth		Chest support pad with lateral guide can be adjusted in height, width and depth, additionally angle adjustable		Back strap – for chest support with lateral guide	
								
	Inside width	Height		Chest width		Chest width		
Size 1	13 cm	27 cm	Size 1	16 – 24 cm	Size 1	16 – 24 cm		
Size 2	16 cm	36 cm	Size 2	18 – 26 cm	Size 2	18 – 26 cm	Size 1	17 cm
Size 3	19 cm	43 cm	Size 3	24 – 34 cm	Size 3	24 – 34 cm	Size 2	28 cm
			Size 4	30 – 42 cm	Size 4	30 – 42 cm	Size 3	39 cm
Price			Upcharge / piece		Upcharge / piece		Price	
☐ Size 1 (4901001)	incl.		☐ Size1 (4901002)	£ 193	☐ Size1 (4901003)	£ 193	☐ Size 1 (4901024)	£ 81
☐ Size 2 (4902001)	incl.		☐ Size2 (4902002)	£ 193	☐ Size 2 (4902003)	£ 193	☐ Size 2 (4902024)	£ 81
☐ Size 3 (4903001)	incl.		☐ Size3 (4903002)	£ 193	☐ Size 3 (4903003)	£ 193	☐ Size 3 (4903024)	£ 81
			☐ Size4 (4904002)	£ 193	☐ Size4 (4904003)	£ 193		

Configuration possibilities - Bracket for buttock and pelvic positioning

Configuration 1:



Bracket for rear support pad (standard)

depth adjustable via quick adjustment

	Depth adjustment Pelvic pad up to buttock pad	Standard configuration
Size 1	11 - 24 cm	for size 1+2
Size 2	15 - 28,5 cm	for size 3+4
		Price
☐ Size 1 (4901005)		incl.
☐ Size 2 (4903005)		incl.

Configuration 2:




Bracket for rear support pad (flip away)

for step-by-step adjustment in the depth via hand crank and quick adjustment, fold-down mechanism for easier transfer (approx. 20 ° (10 cm) downwards)

	Depth adjustment Pelvic pad up to buttock pad	Depth adjustment via Crank	Standard configuration
Size 1	11 - 20 cm	6 cm	for size 1+2
Size 2	15 - 27 cm	6 cm	for size 3+4
			Upcharge / set
☐ Size 1 (4901004)			£ 208
☐ Size 2 (4903004)			£ 208

Configuration 3:



Pelvic frame (with swivel arm)

continuously adjustable in depth and width (depending on the configuration of the pads on page 130), additional depth adjustment by hand crank for gradual straightening of the user, incl. holder for buttocks support pad and back support pad

new


available from
spring 2021.

	Depth adjustment Pelvic pad up to buttock pad	Depth adjustment via Crank	
	11 - 24 cm	6 cm	
			Upcharge / set
☐ for size 1 (4901029)			£ 310
☐ for size 2 (4902029)			£ 310

Configuration options - pelvis positioning

Option 1:

Pelvic support pad with lateral guide + rear support pad



Pelvic support pad with lateral guide width and depth adjustable



	Pelvic width	Depth adjustment
Size 1	16 – 24 cm	5 cm
Size 2	18 – 26 cm	5 cm
Size 3	24 – 34 cm	5 cm
Size 4	30 – 42 cm	5 cm
		Price
☐ Size 1 (49 01 006)		£ 214
☐ Size 2 (49 02 106)		£ 214
☐ Size 3 (49 02 006)		£ 214
☐ Size 4 (49 03 006)		£ 214


Rear support pad



	Width x Height	
Size 1	16x10 cm	
Size 2	20x10 cm	
Size 3	25x14 cm	
		Price
☐ Size 1 (49 01 007)		incl.
☐ Size 2 (49 02 007)		incl.
☐ Size 3 (49 03 007)		incl.

Option 2:

Pelvic positioner + pelvic support pad with lateral guide



Pelvic positioner – depth adjustable



	Width x Height	Depth adjustment
Size 1	16 x 10 cm	5 cm
Size 2	20 x 10 cm	5 cm
Size 3	25 x 14 cm	5 cm
		Price
☐ Size 1 (49 01 008)		£ 83
☐ Size 2 (49 02 008)		£ 83
☐ Size 3 (49 03 008)		£ 83


Pelvic support pad with lateral guide width and depth adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Upcharge / piece	
☐ Size 1 (49 01 010)	£ 133
☐ Size 2 (49 02 110)	£ 133
☐ Size 3 (49 02 010)	£ 133
☐ Size 4 (49 03 010)	£ 133

Option 3:

Pelvic positioner with independent lateral support pads + rear support pad



Pelvic positioner depth adjustable



	Width x Height	Depth adjustment
Size 1	16 x 10 cm	5 cm
Size 2	20 x 10 cm	5 cm
Size 3	25 x 14 cm	5 cm
		Price
☐ Size 1 (49 01 008)		£ 81
☐ Size 2 (49 02 008)		£ 81
☐ Size 3 (49 03 008)		£ 81

Combined spine and pelvic support pads – width, depth and height adjustable



	Pelvic width	Depth adjustment
Size 1	16 – 26 cm	5 cm
Size 2	22 – 32 cm	5 cm
Size 3	32 – 42 cm	5 cm
		Price / pair
☐ Size 1 (49 01 009)		£ 214
☐ Size 2 (49 02 009)		£ 214
☐ Size 3 (49 03 009)		£ 214

Rear support pad



	Width x Height
Size 1	16x10 cm
Size 2	20x10 cm
Size 3	25x14 cm
Price	
☐ Size 1 (4901 007)	incl.
☐ Size 2 (4902 007)	incl.
☐ Size 3 (4903 007)	incl.

Accessories



Therapy tray – made of wood, removable, with edge, incl. height adjustable attachment for chest support pad

	Width	Depth
for size 3+4	70 cm	58 cm
Price		
☐ for size 3+4 (4903 028)		
£ 208		



Therapy tray with edge – made of Perspex, removable, not usable in combination with basket

	Width	Depth
for size 1+2	52 cm	40 cm
for size 3+4	62 cm	50 cm
Price		
☐ for size 1+2 (4901013)		
☐ for size 3+4 (4903 027)		
£ 208		

Accessories

	Basket – removable, not usable in combination with therapy tray	
		Price / pair
☐ for size 1+2 (4900014)		£ 91

	Mounting position of the brake – the brake lever can be mounted under the tray on the left or right, as desired. (standard configuration = under the tray on the right)	
		Price
☐ under the tray on the left (4901015)		incl.
☐ under the tray on the right (4902015)		incl.

Drive wheels

Drive wheels – with drum brake, push rims and black PU tyres		Drive wheels with closed hand rim with drum brake and PU tyres	
			
		for 24" and 28" = black PU tyres for 30" – 40" = grey PU tyres	
		Price / pair	
☐ 24" for size 1 (4900016)	£ 229	☐ 24" for size 1 (4900017)	£ 377
☐ 28" for size 1 (4901016)	£ 229	☐ 28" for size 1 (4901017)	£ 377
☐ 30" for Size 2 (4902016)	£ 279	☐ 30" for size 2 (4902017)	£ 419
☐ 32" for size 2+3 (4903016)	£ 279	☐ 32" for size 2+3 (4903017)	£ 419
☐ 36" for Size 2-4 (4904016)	£ 279	☐ 36" for size 2-4 (4904017)	£ 419
☐ 40" for size 4 (4905016)	£ 329	☐ 40" for size 4 (4905017)	£ 487
☐ 42" for size 4 (4906016)	£ 329		

Drive wheels

	Spoke guard – transparent		
		Price / pair	
	☐ 24" for size 1	(4900018)	£ 120
	☐ 28" for size 1	(4901018)	£ 120
	☐ 30" for Size 2	(4902018)	£ 120
	☐ 32" for size 2+3	(4903018)	£ 175
	Wheel camber adjustment to set the hand rim width		
		Price	
	☐ 0° for size 3+4	(4900019)	incl.
	☐ 3° for size 2-4	(4901019)	incl.
	☐ 6° for size 2-4	(4902019)	incl.
	☐ 9° for size 1-4	(4903019)	incl.
	☐ 12° for size 1+2	(4904019)	incl.
	☐ 15° for size 1	(4905019)	incl.

Accessories

Footplate – height and depth adjustable		Split footplate – individual height-, depth-, width- and angle adjustable already included in basic model „with abduction“	
			
		Heel edges with ankle straps – angle, width and depth adjustable, for attachment on the footplate/ split footplate with magnetic lock.	
			
		Front width	
		Size 1	8.5 – 11 cm
		Size 2	10 – 12.5 cm
		Price / pair	
☐ for size 1 (4901020)	incl.	☐ for size 1 (4901021)	£ 221
☐ for size 2 (4902020)	incl.	☐ for size 2 (4902021)	£ 221
☐ for size 3 (4903020)	incl.	☐ for size 3+4 (4903021)	£ 221
☐ for size 4 (4904020)	incl.	☐ Size 1 (4901027)	£ 65
		☐ Size 2 (4902027)	£ 65

☐ **Order**

☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____

Customer:

Notes:

Reference:

The standing frame – from prone to supine.


- Basic frame with height and depth adjustable footplate
- Angular adjustment of the central column via gas pressure spring from (-) 15° – 30°
- Incl. chest support pad with lateral guide and rear support pad
- Height and depth adjustable pelvis frame, pelvis position adjustable via a crank
- Incl. three-dimensionally adjustable knee support pads with patella cutout
- Incl. therapy tray
- 4 castors with locking devices (75 mm) and total lock

new Frame colour
Turquoise blue



Please send us the following measurements for optimum presetting of the standing frame

Settings



The following measurements should be pre-set on the standing frame:

A Body size _____

B Sternum height _____

C Elbow height _____

D Crotch height _____

E Chest width _____

F Pelvis width _____

G Knee width _____


H Knee height _____

☐ **48 00 020**

Price

incl.

Dimensions - basic model

		Size 1	Size 2	Size 3	
	Body size	75 - 110 cm	100 - 140 cm	120 - 160 cm	
A	Height from footplate to tray	48.5 - 80.5 cm	71 - 95 cm	80 - 112 cm	
B	Pelvis depth	12 - 22 cm	16 - 28 cm	20 - 35 cm	
	Depth adjustment of pelvis frame (via crank)	6 cm	6 cm	6 cm	
C	Total length	84 cm	90 cm	100 cm	
D	Total width	57 cm	65 cm	71 cm	
E	Abduction split footplate	0° - 27°	0° - 27°	0° - 27°	
F	Pelvis height	45 - 60 cm	57 - 80 cm	65 - 90 cm	
G	Sternum height	55 - 85 cm	75 - 110 cm	95 - 125 cm	
	Max. load	40 kg	60 kg	80 kg	

Basic model

tim. size 1 for users with a body size 75 - 110 cm		tim. size 2 for users with a body size 100 - 140 cm		tim. size 3 for users with a body size 120 - 160 cm	
Price		Price		Price	
☐ 48 01 100	£ 1,650	☐ 48 02 100	£ 1,650	☐ 48 03 100	£ 1,750

Frame colours

Colour selection for size 1		Colour selection for size 2		Colour selection for size 3		Special colour (on request)
						Please provide the requested colour here: RAL _____
Price		Price		Price		Price
☐ Turquoise blue (92 00 016)	incl.	☐ Turquoise blue (92 00 016)	incl.	☐ Turquoise blue (92 00 016)	incl.	☐ 48 00 090
☐ Pure white (92 00 100)	incl.	☐ Pure white (92 00 100)	incl.	☐ Pure white (92 00 100)	incl.	
☐ Melon yellow (92 00 006)	incl.	☐ Dormant Red (92 00 003)	incl.	☐ Lazurite blue metallic (92 00 010)	incl.	

Coated covers

Flexi Black		Flexi Light grey		Fabric choice	
				Fabric choice – please enter here the desired fabric out of our fabric collection: Art.-No. (see page 156) _____	
Price		Price		Price	
☐ 91 00 043	incl.	☐ 91 00 050	incl.	☐ 91 00 999	incl.


Lateral pads

Chest support pad with lateral guide can be adjusted in height, width and depth		Chest support pad with lateral guide can be adjusted in height, width and depth, additionally angle adjustable		Chin pad – height adjustable		Back strap	
							
Chest width		Chest width		Inside width		Pad width	Pad height
Size 1	16 - 24 cm	Size 1	16 - 24 cm	Size 1	6 cm	Size 1	17 cm
Size 2	18 - 26 cm	Size 2	18 - 26 cm	Size 2	8 cm	Size 2	28 cm
Size 3	24 - 34 cm	Size 3	24 - 34 cm	Size 3	10 cm	Size 3	39 cm
Size 4	30 - 42 cm	Size 4	30 - 42 cm				
Price		Upcharge / piece		Price		Price	
☐ Size 1 (48 01 101)	incl.	☐ Size 1 (48 01 102)	£ 68	☐ Size 1 (48 01 012)	£ 128	☐ Size 1 (48 01 013)	£ 77
☐ Size 2 (48 02 201)	incl.	☐ Size 2 (48 02 202)	£ 68	☐ Size 2 (48 02 012)	£ 128	☐ Size 2 (48 02 013)	£ 77
☐ Size 3 (48 02 101)	incl.	☐ Size 3 (48 02 102)	£ 68	☐ Size 3 (48 03 012)	£ 128	☐ Size 3 (48 03 013)	£ 77
☐ Size 4 (48 03 101)	incl.	☐ Size 4 (48 03 102)	£ 68				

Combinations - pelvic frame

Option 1:

Pelvis support pad with lateral guide (mounted at the front of the pelvic frame) + rear support pad (mounted at the back of the pelvic frame)



Pelvis support pad with lateral guide width adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Price	
☐ Size 1 (48 01 005)	£ 198
☐ Size 2 (48 02 105)	£ 198
☐ Size 3 (48 02 005)	£ 198
☐ Size 4 (48 03 005)	£ 198


Rear support pad



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price	
☐ Size 1 (48 01 003)	incl.
☐ Size 2 (48 02 003)	incl.
☐ Size 3 (48 03 003)	incl.

Option 2:

Pelvic positioner (mounted at the front of the frame) + pelvic support pad with lateral guide (mounted at the back of the frame)



Pelvic positioner



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price	
☐ Size 1 (48 01 004)	£ 88
☐ Size 2 (48 02 004)	£ 88
☐ Size 3 (48 03 004)	£ 88


Pelvic support pad with lateral guide width adjustable



	Pelvic width
Size 1	16 – 24 cm
Size 2	18 – 26 cm
Size 3	24 – 34 cm
Size 4	30 – 42 cm
Upcharge / piece	
☐ Size 1 (48 01 006)	£ 198
☐ Size 2 (48 02 106)	£ 198
☐ Size 3 (48 02 006)	£ 198
☐ Size 4 (48 03 006)	£ 198

Option 3:

Pelvic positioner with independent lateral support pads (mounted at the front of the frame) + rear support pad (mounted at the rear of the frame)



Pelvic positioner



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price	
☐ Size 1 (48 01 004)	£ 88
☐ Size 2 (48 02 004)	£ 88
☐ Size 3 (48 03 004)	£ 88

Lateral support pads – can be independently adjusted in depth, height and width



	Pelvic width
for size 1	16 – 26 cm
for size 2	22 – 32 cm
for size 3	32 – 42 cm
Price / pair	
☐ for size 1 (48 01 007)	£ 198
☐ for size 2 (48 02 007)	£ 198
☐ for size 3 (48 03 007)	£ 198

Rear support pad



	Width x Height
Size 1	16 x 10 cm
Size 2	20 x 10 cm
Size 3	25 x 14 cm
Price	
☐ Size 1 (48 01 003)	incl.
☐ Size 2 (48 02 003)	incl.
☐ Size 3 (48 03 003)	incl.


Support pads				Accessories			
Knee support pads – individual three-dimensional adjustment, with patella recess		Knee support pads – individual three-dimensional adjustment, with adjustable attachment brackets (knee flexion)		Back support pad – mounted at the back of the pelvis frame, height and angle adjustable		Headrest – height, depth and angle adjustable, only for use together with the back support pad	
	Inside width		Inside width		Width	Height	
Size 1	8 cm	Size 1	8 cm	for size 1	20 cm	23 cm	
Size 2	10 cm	Size 2	10 cm	for size 2	20 cm	23 cm	
Size 3	12 cm	Size 3	12 cm	for size 3	25 cm	28 cm	
Price / pair		Upcharge / pair		Price		Price	
☐ Size 1 (48 01 008)	incl.	☐ Size 1 (48 01 009)	£ 68	☐ for size 1 (48 01 010)	£ 218	☐ for size 1 (48 01 111)	£ 175
☐ Size 2 (48 02 008)	incl.	☐ Size 2 (48 02 009)	£ 68	☐ for size 2 (48 02 010)	£ 218	☐ for size 2 (48 02 111)	£ 175
☐ Size 3 (48 03 008)	incl.	☐ Size 3 (48 03 009)	£ 68	☐ for size 3 (48 03 010)	£ 218	☐ for size 3 (48 03 111)	£ 175

Accessories			
Therapy tray	Therapy tray – with bowl (incl. cover)	Therapy tray – with adjustable armrests	Therapy tray – with adjustable armrests and bowl (incl. cover)
			
Price	Upcharge / piece	Upcharge / piece	Upcharge / piece
☐ for size 1 (48 01 114)	incl.	☐ for size 1 (48 01 115)	£ 128
☐ for size 2 (48 02 114)	incl.	☐ for size 2 (48 02 115)	£ 128
☐ for size 3 (48 03 114)	incl.	☐ for size 3 (48 03 115)	£ 128
		☐ for size 1 (48 01 116)	£ 113
		☐ for size 2 (48 02 116)	£ 113
		☐ for size 3 (48 03 116)	£ 113
		☐ for size 1 (48 01 117)	£ 244
		☐ for size 2 (48 02 117)	£ 244
		☐ for size 3 (48 03 117)	£ 244

Accessories			
Padded tray cover without armrests covered with black Flexi cover, height: 1 cm	Pad for therapy tray with armrests covered with black Flexi cover, height: 1 cm	Angular adjustment of the center column by gas spring from (-) 15° - 30° – Standard configuration = with foot pedal, if nothing else stated	
			
Price	Price	Price	
☐ for size 1 (48 01 018)	£ 92	☐ for size 1 (48 01 019)	£ 113
☐ for size 2 (48 02 018)	£ 92	☐ for size 2 (48 02 019)	£ 113
☐ for size 3 (48 03 018)	£ 92	☐ for size 3 (48 03 019)	£ 113
		☐ by foot pedal on the base frame (48 00 025)	incl.
		☐ by manual release under the tray (48 01 025)	incl.

Accessories			
Footplate – height and depth adjustable	Split footplate – height, width and angle adjustable (abduction of legs)	Heel edges with ankle straps – angle, width and depth adjustable, for attachment on the footplate/ split footplate with magnetic lock	
			
		Front width	
		Size 1	8.5 – 11 cm
		Size 2	10 – 12.5 cm
Price	Upcharge / pair	Price / pair	
☐ for size 1 (48 01 120)	incl.	☐ Size 1 (48 01 024)	£ 66
☐ for size 2 (48 02 120)	incl.	☐ Size 2 (48 02 024)	£ 66
☐ for size 3 (48 03 120)	incl.		

Body measurements - standing



The following measurements can be preset:

A Body size _____

E Chest width _____

B Sternum height _____

F Pelvis width _____

C Elbow height _____

G Knee width _____

D Crotch height _____

H Knee height _____



malte.

The compact and versatile posterior walker.

Further information can be found from page 30.



Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____


Reference:
The comfortable toilet and shower chair.

- Anatomically-shaped seat pan made of PE, colour steel blue
- Subframe made of stainless steel
- Tilt-in-space adjustment via gas pressure spring:
130N (size 0), 200N (size 1-3) and 2x 130N (size 4)
- Mechanical height adjustment
- 4 dual castors (100 mm) with locking devices
- Incl. pelvic harness and PU pad


Dimensions - basic model

		Size 0	Size 1	Size 2	Size 3	Size 4
A	Seat width (rear)	26 cm	31 cm	36 cm	41 cm	46 cm
B	Seat depth	31 cm	36 cm	41 cm	46 cm	48 cm
C	Seat shell width (front)	31 cm	36 cm	41 cm	46 cm	51 cm
D	Back height	49 cm	51 cm	56 cm	61 cm	65 cm
E	Tilt-in-space	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°	(-)5° - 30°
F	Can be positioned over a toilet	47 - 57 cm	45 - 55 cm	42 - 52 cm	42 - 52 cm	42 - 52 cm
G	Total width	56 cm	61 cm	63 cm	71.5 cm	74 cm
H	Total length	65 cm	65 cm	65 cm	70 cm	70 cm
I	Max. lower leg length	38 cm	38 cm	48 cm	48 cm	48 cm
J	Front edge of hygiene opening to back edge of subframe	47 cm	49 cm	51 cm	54 cm	57 cm
K	Front opening width	10 cm	10 cm	10 cm	10 cm	10 cm
L	Rear opening width	12 cm	12 cm	13.5 cm	13.5 cm	13.5 cm
M*	Rear opening width	-	-	18.5 cm	18.5 cm	18.5 cm
	Max. load	60 kg	80 kg	100 kg	100 kg	150 kg
	Weight	28 kg	30 kg	33 kg	35 kg	38 kg

* only applies to sizes 2 b, 3 b and 4 b



Basic model				
boris. size 0	boris. size 1	boris. size 2	boris. size 3	boris. size 4
Price	Price	Price	Price	Price
☐ 52 00 100/UK £ 1,685	☐ 52 01 100/UK £ 1,685	☐ Size 2a (52 02 100/UK) £ 1,785 ☐ Size 2b (52 02 200/UK) £ 1,785	☐ Size 3a (52 03 100/UK) £ 1,820 ☐ Size 3b (52 03 200/UK) £ 1,820	☐ Size 4a (52 04 100/UK) £ 1,920 ☐ Size 4b (52 04 200/UK) £ 1,920

Accessories				
PU headrest – blue, adjustable in height and depth with stainless steel bracket 	PU insert pad – two-piece for opening form A 	PU insert pad – two-piece for opening form B 	Tray – with padded arm rest segment and holder 	
Price	Price	Price	Price	Price
☐ 52 00 001 £ 148	☐ for size 0 (52 00 101) incl. ☐ for size 1 (52 01 101) incl. ☐ for size 2 (52 02 101) incl. ☐ for size 3 (52 03 101) incl. ☐ for size 4 (52 04 101) incl.	☐ for size 2 (52 02 201) incl. ☐ for size 3 (52 03 201) incl. ☐ for size 4 (52 04 201) incl.	☐ for size 0 (52 00 012) £ 204 ☐ for size 1 (52 01 002) £ 204 ☐ for size 2 (52 02 002) £ 204 ☐ for size 3 (52 03 002) £ 204 ☐ for size 4 (52 04 002) £ 204	

Accessories																								
Toilet pan – incl. lid and bracket as assembly set	Abduction wedge – integrated splash guard (Note: a subsequent installation with size 4 must be done in the factory!)		Footplate – height, knee and foot angle adjustment, also folding		PU-Calf plate – height adjustable (only for use in combination with the footplate)																			
																								
		<table><tr><td></td><td>Width</td><td>Height from seat</td></tr><tr><td>for size 0+1</td><td>7 cm</td><td>14 cm</td></tr><tr><td>for size 2-4</td><td>8 cm</td><td>17.5 cm</td></tr></table>		Width	Height from seat	for size 0+1	7 cm	14 cm	for size 2-4	8 cm	17.5 cm		<table><tr><td></td><td>Adjustable lower leg length</td></tr><tr><td>for size 0+1</td><td>30 – 38 cm</td></tr><tr><td>for size 2-4</td><td>35 – 48 cm</td></tr></table>		Adjustable lower leg length	for size 0+1	30 – 38 cm	for size 2-4	35 – 48 cm	<table><tr><td>Height</td><td>15 cm</td></tr><tr><td>Width</td><td>38 cm</td></tr></table>	Height	15 cm	Width	38 cm
	Width	Height from seat																						
for size 0+1	7 cm	14 cm																						
for size 2-4	8 cm	17.5 cm																						
	Adjustable lower leg length																							
for size 0+1	30 – 38 cm																							
for size 2-4	35 – 48 cm																							
Height	15 cm																							
Width	38 cm																							
Price / set		Price		Price																				
☐ for size 0+1+3 (52 00 102)	£ 95	☐ for size 0+1 (52 01 103/UK)	£ 99	☐ for size 0+1 (52 01 004/UK)	£ 184																			
☐ for size 2+4 (52 02 102)	£ 95	☐ for size 2+4 (52 02 103/UK)	£ 99	☐ for size 2-4 (52 02 004/UK)	£ 184																			
		☐ for size 3 (52 03 103/UK)	£ 99																					
				☐ 52 01 009	£ 115																			

Accessories		Harness holder			
Fixation harness – made of mesh fabric incl. strap holder for shoulder straps, colour: blue 	Pelvic harness – made of mesh fabric, colour: blue (standard configuration = pelvic harness size 1 for boris. size 0+1, pelvic harness size 2 for boris. size 2 and pelvic harness size 3 for boris. size 3+4, if nothing else stated) 	Strap holder – for shoulder straps when using other fixation aids apart from the adjacent positioning jacket 			
Price / set	Price	Art.-No.	Use	Max. shoulder width	Price
☐ Size 1 (52 01 010/UK) £ 129 ☐ Size 2 (52 02 010/UK) £ 129 ☐ Size 3 (52 03 010/UK) £ 129	☐ Size 1 (52 01 011) incl. ☐ Size 2 (52 02 011) incl. ☐ Size 3 (52 03 011) incl.	☐ 52 01 005 ☐ 52 01 008	with headrest without headrest	30 cm 30 cm	£ 62 £ 62

Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____

Signature: _____



Reference:
The versatile bathing and shower cradle.

- Black, powder-coated frame with removable, washable cover
- Incl. back angle, knee angle and head angle adjustment
- Incl. chest strap


Dimensions - basic model

		Size 1	Size 2	Size 3	Size 4
A	Seat width	30 cm	35 cm	35 cm	40 cm
B	Seat depth	30 cm	34 cm	38 cm	44 cm
C	Back height	45 cm	61 cm	71 cm	82 cm
D	Standing space length	69 cm	74 cm	88 cm	93 cm
E	Lower leg length	19 cm	20 cm	28 cm	28 cm
F	Total width	36 cm	41 cm	41 cm	46 cm
	Back angle	90° - 180°	90° - 180°	90° - 180°	90° - 180°
	Knee angle	0° - 90°	0° - 90°	0° - 90°	0° - 90°
	Headrest angle	0° - 90°	0° - 90°	0° - 90°	0° - 90°
	Seat height	0° / 5 cm	0° / 5 cm	0° / 5 cm	0° / 5 cm
	Max. load	35 kg	45 kg	55 kg	75 kg
	Weight	4 kg	4.1 kg	4.3 kg	4.8 kg

* in combination with the cover with pelvic mould depression



Basic model

basti. size 1		basti. size 2		basti. size 3		basti. size 4 Back angle adjustment via stainless steel ratchet joint	
Price		Price		Price		Price	
☐ 56 01 000	£ 515	☐ 56 02 000	£ 565	☐ 56 03 000	£ 615	☐ 56 04 000	£ 715

Accessories

Cover with hollow for pelvis and adjustable tube edge




Head guide pad – adjustable with Velcro



Upcharge / piece				Price	
☐ for size 1 (56 01 001)		£ 88	☐ for size 1 (56 01 002)	£ 88	£ 105
☐ for size 2 (56 02 001)		£ 88	☐ for size 2 (56 02 002)	£ 88	£ 105
☐ for size 3 (56 03 001)		£ 88	☐ for size 3 (56 03 002)	£ 88	£ 105
☐ for size 4 (56 04 001)		£ 88	☐ for size 4 (56 04 002)	£ 88	£ 105

Accessories

Additional retainer strap – adjustable with Velcro for the chest or pelvis





basti. shower cradle subframe – mobile with 4 castors with brakes



Price				Price	
☐ for size 1 (56 01 003)		£ 88	☐ 56 00 200		£ 515
☐ for size 2+3 (56 03 003)		£ 88			
☐ for size 4 (56 04 003)		£ 88			

Dimensions - subframe

A	Length	92 cm
B	Width	57 cm
C	Seat height with subframe	61" / 66 cm
	* in combination with the cover with pelvic mould depression	



Customer:
Notes:
☐ **Order**
☐ **Cost estimate**

Customer No.: _____

Contact: _____

Date: _____


Signature: _____

Reference:
The lightweight and compact bathing seat.

- Aluminium frame with removable cover
- Back section with adjustable angle
- Incl. pelvic harness


Dimensions - basic model

		bibi.	bibi. with long back
A	Seat depth	30 cm	30 cm
B	Seat width	36 cm	36 cm
C	Back height	45 cm	55 cm
D	Total length	60 cm	60 cm
	Back angle	90° - 180°	90° - 180°
	Max. load	30 kg	30 kg
	Weight	2.3 kg	2.3 kg


Basic model


bibi.		bibi. with long back	
	Price		Price
☐ 55 01 000	£ 280	☐ 55 01 001	£ 280

Accessories
Additional pelvic harness


☐ 55 01 002	£ 80
---	------

Body measurements - hygiene

63



The following measurements can be preset:

A Body size	_____	F Knee width	_____
B Sternum height	_____	G Knee height	_____
C Inside leg length	_____	H Seat depth	_____
D Chest width	_____	I Elbow height	_____
E Pelvis width	_____	J Back height	_____

Cleaning + disinfection.

Introduction to using our products.

Please note that **people with special needs, reading difficulties and impaired vision** must ask a parent, supervisor or technical advisor to read out or explain this technical documentation so that safe usage of our products can be guaranteed.

We would like to point out that only **spare parts and accessories** from our company may be used and installed so as not to endanger the safety of the users. In individual cases, the combination of our products with those of other manufacturers may require clarification with your specialist retailer or with us as the manufacturer.

Settings and adjustments to the product or accessories may only be made by people who have been given the necessary instructions by a medical product advisor. Please ensure that none of the user's extremities are in the respective area when making adjustments of any kind to minimise the risk of injury

The respective product should be reduced to its smallest dimensions **for transportation**. For example, reduce the seat, the armrests, the back etc. to the smallest dimensions. Remove the parts of the product that can be easily removed, e.g. the wheels of the **todd**. standing frame. Fold up the foot support on products with this feature. Please see the User Manual belonging to the respective product for more detailed information.

Instructions for use on the respective products are available on www.schuchmann.co.uk/media-library.

If applicable, the **product-specific cleaning guidelines** in the respective instructions for use are to be taken into consideration.

Cleaning guidelines for the different materials.

The care and cleaning of our products should be carried out at regular intervals to ensure that they remain functional over their useful life. The choice of suitable cleaning agents is crucial.

Wood: Wooden parts should be cleaned with a damp cloth or similar. When using detergents, please consult the manufacturer's instructions regarding detergent concentration. Caustic household detergents or corrosive substances are not suitable. Never leave natural wood in a damp or wet condition over long periods, as it will take on water and can change in shape (swelling). Therefore, only wipe and clean damp (i.e. with wrung-out cloths). Test a small, non-visible area for compatibility.

Powder-coated parts: Aluminium parts should be cleaned with a sponge or damp cloth using normal antibacterial household detergents.

Stainless steel: Stainless steel parts should be cleaned using a sponge or damp cloth with a mild soap. Bleach must not be used.

Covers / Fabric: Covers and fabric should be washed at 40° C in the washing machine using normal detergent. Do not tumble dry.

Plastics: Plastics should be cleaned using a sponge or damp cloth with mild soap. Bleach must not be used.

PU (polyurethane): PU parts should only be cleaned using a sponge or damp cloth with a mild soap. Bleach must not be used.

Should you have any
queries, you are welcome
to contact us directly!



